

# IC-Flow™ Imaging System

Diagnostic  
**Green**



IC-Flow™ Imaging System is the world's first handheld fluorescence imaging device with integrated display of the fluorescence image.



# IC-Flow procedure


## Fluorescence principle

The fluorescent dye is dissolved and injected. It emits a fluorescence light when it is excited by the near-infrared (NIR) light source of the IC-Flow Camera.

The resulting fluorescent light is detected by the IC-Flow Camera and displayed on the IC-Flow Monitor or internal touchscreen.

Fluorescent images can be obtained at a maximum 1.5 cm below the tissue surface.

The IC-Flow takes advantage of the NIR content of ambient lighting so that the body contours are visible. This is important for image interpretation during a procedure.





## Unique user benefits

- Flexible configurations to meet user-specific needs
  - Integrated Operating Room System
  - Portable Point of Care System
- Handheld and easy to use
- Benefits of LED compared with laser light source
  - No maintenance
  - Minimal service requirements
- Operated via camera or touchscreen
- Fast boot up time
- Backed by Diagnostic Green
  - Global supplier of ICG
  - 20 years experience in Perfusion Imaging
  - IC-View camera -> PDE camera -> IC-Flow

## The use of IC-Flow Imaging System

1. After being switched on, the system is ready to be used within seconds
2. Using the Green Balance ICG Reference Card -> functional check and camera settings adjustments
3. Fluorescent dye is dissolved and injected
4. IC-Flow Camera is held at a distance of 15 - 20 cm to the tissue
5. Display of fluorescence image on monitor and/or touchscreen
6. Recording a video or taking a snapshot
7. Transfer of stored data to an USB stick (demonstration purpose)



# IC-Flow Configurations


## Point of Care System

- Compact and easily portable
- Fluorescence image on touchscreen
- No additional equipment required
- Data transfer to USB stick (demo purpose)



## Operating Room System

- Cart system for organization
- Integrated system with monitor
- Fluorescence image on touchscreen & monitor
- External recorder for official recording with patient data



# IC-Flow & Accessories



**IC-Flow™ Monitor**  
(Art. No. PC6224)



**IC-Flow™ Imaging System**  
(Art. No. PC6200)




**IC-Flow™ Video Recorder**  
(Art. No. PC225)



**IC-Flow™ Cart**  
(Art. No. PC220)



**Green Balance ICG Reference Card**  
(Art. No. PV6210U)



**IC-Flow™ Camera Cover**  
(Art. No. PV6240S)

IC-Flow™  
Monitor

IC-Flow™  
Cart

# Benefits - IC-Flow method

## Why IC-Flow method?

Using the IC-Flow allows a **real-time visualization** of the injected dye in the tissue.

Due to the fast preparation time and the short period of use the method can be performed **before, while and after the treatment procedure.**

The usage of this imaging method can be **learned easily.**

In comparison with some other imaging methods the IC-Flow procedure causes no **harmful radiation.**





SLN detection procedure with IC-



# Indications

## Plastic reconstruction surgery


- Perfusion visualization of flaps (patency/control check)
- Visualization of intraoperative anastomosis



*Dr. Ulf Dornseifer, Klinikum Bogenhausen, Munich Germany.*

## Secondary Lymphedema


- Evaluation support of lymphatic function
- Assessment of lymphatic flow velocity




*JP. Belgrado, L. Vandermeeren, Université libre de Bruxelles, Brussels, Belgium.*

## Oncology (Gynecology/ Dermatology)

- Visualization of the lymphatic vessels
- Detection of Sentinel Lymph Nodes (breast/skin cancer)



*Dr. Bianca Baican, AGAPLESION Markus Krankenhaus, Frankfurt, Germany.*



h IC-Flow, Universitätsklinikum Tübingen, Germany



**Manufacturer & Distributor**

**Diagnostic Green GmbH**  
Otto-Hahn-Straße 20  
85609 Aschheim-Dornach  
Germany

Telephone: +49 (0) 89.1241.477.20  
Fax: +49 (0) 89.1241.477.29  
info@diagnosticgreen.com  
www.diagnosticgreen.com

