INSTRUCTIONS FOR USE FOR:

T I S S U E
R E I N F O R C E M E N T

en

English

hu

Magyar

cs

Čeština

it

Italiano

da

Dansk

no

Norsk

nl

Nederlands

pl

Polski

fi

Suomi

pt

Português

fr

Français

es

Español

de

Deutsch

sv

Svenska

el

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tr

Türkçe

INSTRUCTIONS FOR USE GORE® BIO-A® TISSUE REINFORCEMENT

Carefully read all instructions prior to use. Observe all instructions, warnings, and precautions noted throughout. Failure to do so may result in complications.

PRODUCT DESCRIPTION

As packaged, the GORE® BIO-A® Tissue Reinforcement is a tailorable, bioabsorbable material comprised of synthetic bioabsorbable poly (glycolide:trimethylene carbonate) copolymer (PGA:TMC). The GORE® BIO-A® Tissue Reinforcement is used to reinforce soft tissue during the phases of wound healing by filling soft-tissue deficits.

The GORE® BIO-A® Tissue Reinforcement is designed with a textured porous fibrous web surface on both surfaces. The device elicits a physiological response, which fills the deficit with native tissue and gradually absorbs the device. It is recommended that the device be placed next to well-vascularized tissue where ingrowth is desired. Degraded via a combination of hydrolytic and enzymatic pathways, the PGA:TMC copolymer has been found to be both biocompatible and non-immunogenic. In vivo studies with this copolymer indicate the bioabsorption process should be complete by six to seven months¹.

INDICATIONS FOR USE

The GORE® BIO-A® Tissue Reinforcement is intended for use in the reinforcement of soft tissue. Examples where the GORE® BIO-A® Tissue Reinforcement may be used are hiatal and ventral hernia repair as suture-line reinforcement.

CONTRAINDICATIONS

The GORE® BIO-A® Tissue Reinforcement is contraindicated for use in reconstruction of cardiovascular defects.

Because GORE® BIO-A® Tissue Reinforcement is absorbable, it is contraindicated for use in patients requiring permanent support from the device.

WARNINGS AND PRECAUTIONS

Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient (SEE ADVERSE EVENTS).

WARNINGS

- Do not reuse. Do not use if the package is open or damaged, past the expiration date, or if it is suspected that the sterility of the device has been compromised. This may lead to infection or affect product performance and related serious potential harms could result. Due to the absorbable nature of the device these conditions may result in early device degradation.
- GORE® BIO-A® Tissue Reinforcement placement with unnecessary contact with small and/or large bowel should be avoided because it may result in adhesions. Clinical data on abdominal surgery, with or without surgical mesh, has shown that adhesion formation may occur and may result in complications including pain, bowel obstruction, and additional intervention including surgery.
- Failure to adequately place or fixate the device, or placing the device with excessive folds may result in inadequate tissue ingrowth or recurrence of the defect. This may lead to bowel obstruction, dyspareunia, dysphagia, GERD recurrence, erosion or extrusion and related harms, fistula, infection, pain, secondary procedure, seroma or hematoma and related harms, tissue ischemia, urinary problems, vaginal bleeding or discharge, vaginal dysfunction and wound dehiscence.
- Use of this product in applications other than those indicated, including but not limited to those listed below, has the potential for serious complications.

- The safety and effectiveness of this device for use in transvaginal insertion techniques has not been established. Transvaginal insertion techniques, which expose the surgical mesh to the vaginal flora, can increase the risk of contamination leading to bacterial colonization of the surgical mesh. Prolonged exposure to bacteria may necessitate surgical mesh removal. If considering use of surgical mesh for transvaginal stress urinary incontinence and / or pelvic organ prolapse procedures, consult with clinical guidelines for recommendations (e.g. surgical mesh material including device surface area, construction, pore size).
- The safety and effectiveness of this device has not been established for use in transvaginal pelvic organ prolapse repair. There is insufficient evidence to assure that the probable benefits of surgical mesh used in transvaginal pelvic organ prolapse repair outweigh its probable risks, therefore this use is not recommended for this device.

Use of this product in transvaginal or transabdominal pelvic organ prolapse and stress urinary incontinence is not within the indication. Further, refer to clinical guidelines that have recommended that use of any surgical mesh in pelvic organ prolapse and stress urinary incontinence procedures could result in device-related or procedure-related complications, including, but not limited to, bowel obstruction and related harms, dyspareunia, fistula, hematoma, infection, irritation or inflammation, mesh contraction, mesh erosion, mesh exposure, organ injury, pelvic pain, recurrence, urinary problems, vaginal dysfunction, or

vaginal bleeding or discharge, which may require additional surgery including device removal. Complete removal of mesh may not be possible and may not result in complete resolution of complications.

- The safety and effectiveness of this device has not been established for use in patients who have a surgical site infection that cannot be treated successfully prior to device placement (e.g. use of surgical mesh in an infected surgical site). Use of any surgical mesh in the presence of infection could result in potentially serious patient harms and additional intervention, including surgery.

PRECAUTIONS

- The GORE® BIO-A® Tissue Reinforcement should only be used by physicians trained in soft tissue reinforcement according to the indicated uses.
- Implanting of multiple devices on top of each other has not been evaluated for use with GORE® BIO-A® Tissue Reinforcement.
- As with any mesh, use of GORE® BIO-A® Tissue Reinforcement in the presence of contamination may result in dyspareunia, fever, infection, irritation or inflammation, pain, and wound dehiscence, and may require removal of the mesh if infection occurs.

ADVERSE EVENTS

POTENTIAL CLINICAL AND DEVICE ADVERSE EVENTS

Possible adverse events with the use of any tissue deficiency prosthesis or soft tissue reinforcement surgical procedures may include, but are not limited to, additional intervention including surgery, adhesions and related harms, bowel obstruction, defect recurrence and related harms, dysphagia, erosion or extrusion and related harms, exposure or protrusion and related harms, fever, fistula, GERD recurrence, infection, irritation or inflammation, pain, paresthesia, seroma or hematoma and related harms, tissue ischemia, and wound dehiscence.

CLINICAL AND DEVICE ADVERSE EVENT REPORTING

Any adverse event involving the GORE® BIO-A® Tissue Reinforcement should be reported to the manufacturer and the country specific regulatory authorities immediately. To report an event to W. L. Gore & Associates, email: medcomplaints@wlgore.com or contact:

US:

Phone: +1.800.528.1866 or +1.928.864.4922,

Fax: +1.928.864.4364

EMEA:

Phone: +49 89 4612 3440,

Fax: +49 89 4612 43440

To report an adverse event to W. L. Gore & Associates in all other markets, please use the email above (medcomplaints@wlgore.com) or call the US phone number (+1.800.528.1866 or +1.928.864.4922).

HOW SUPPLIED

CONTENTS

The GORE® BIO-A® Tissue Reinforcement is packaged in a sterile inner pouch with a MINIPAX® desiccant pouch. The MINIPAX® desiccant pouch included in the device package is not for implantation. If the MINIPAX® desiccant pouch has been compromised, discard product.

STERILITY

- The GORE® BIO-A® Tissue Reinforcement is provided STERILE and non-pyrogenic, and is sterilized by gamma radiation. Do not resterilize.
- Do not use after the "use by" (expiration) date printed on the label.
- The GORE® BIO-A® Tissue Reinforcement is designed for single use only; do not reuse device. Gore does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination which may lead to patient harms. (See WARNINGS).

STORAGE AND HANDLING

- Store in a dry place. Avoid exposing the package or device to extreme hot or cold temperatures.
- Do not use the GORE® BIO-A® Tissue Reinforcement if the inner or outer pouch is compromised or if the device is damaged.
- To open the package, peel the seal apart from the upper portion of the inner and outer pouch.
- Handle and dispose of the device and packaging, taking into account any infectious or microbial hazard risks, with the necessary precautionary measures in accordance with acceptable medical practice and with applicable local, state, and federal laws and regulations.
- See WARNINGS AND PRECAUTIONS for additional considerations specific to storage and handling.

For additional Storage and Handling information, see HOW SUPPLIED.

DIRECTIONS FOR USE

The GORE® BIO-A® Tissue Reinforcement should only be used by physicians experienced in abdominal wall soft tissue reinforcement.

PATIENT SELECTION AND RECOMMENDED DEVICE SIZING

A physician should determine if a particular patient and anatomy is suitable for proper use of the GORE® BIO-A® Tissue Reinforcement. See WARNINGS AND PRECAUTIONS for additional considerations specific to the procedure.

The device should be properly sized by the physician to fit the patient's anatomy.

DEVICE PREPARATION

- Check information on product labeling for correct configuration, size, and expiration date.
- Open the carton at the flap on the top of the package and remove the pouch.
- Inspect the package and do not use device if the package is damaged or if it is suspected sterility has been compromised.
- Use sterile technique to peel the seal apart from the upper portion of the inner and outer pouch and remove the device.
- Use clean, sterile gloves or atraumatic sterile instruments when handling the GORE® BIO-A® Tissue Reinforcement.
- Using aseptic technique, trim the GORE® BIO-A® Tissue Reinforcement to the desired size using sharp surgical scissors.
- GORE® BIO-A® Tissue Reinforcement does not require pre-wetting. If desired, the device may be wetted with sterile saline or similar standard irrigation solution prior to placement.

IMPLANT PROCEDURE

- Insert the device through an open incision or introduce through a trocar. Ease of trocar introduction may vary depending on rolled device size and trocar availability. Excessive force should be avoided during trocar introduction.
- Place the textured porous fibrous surfaces adjacent to well-vascularized tissue where ingrowth is desired. Ensure the device lays flat without excessive folds. SEE WARNINGS.
- Implanting of multiple devices on top of each other has not been evaluated for use with GORE® BIO-A® Tissue Reinforcement.
- Confirm proper placement of device. If the GORE® BIO-A® Tissue Reinforcement placement is not optimal, reposition device.
- Fixate GORE® BIO-A® Tissue Reinforcement. The device may be sutured to host tissue for stabilization using absorbable or permanent sutures.
- Fixation methods other than those described above have not been evaluated for use with GORE® BIO-A® Tissue Reinforcement.
- The physician should reference the manufacturer's instructions for use and any supporting clinical literature regarding any potential adverse reactions related to fixation devices.

See WARNINGS AND PRECAUTIONS for additional considerations specific to the procedure.

MR SAFETY INFORMATION

The GORE® BIO-A® Tissue Reinforcement is MR Safe.

REFERENCES

¹Katz AR, Mukherjee DP, Kaganov AL, Gordon S. A new synthetic monofilament absorbable suture made from polytrimethylene carbonate. *Surgery, Gynecology & Obstetrics* 1985;161(3):213-222.

DEFINITIONS

[EC REP] Authorised Representative in the European Community

[REF] Catalogue Number

 Caution

R Only CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

 Consult Instructions for Use

 Date of Manufacture

 Do Not Expose to Extreme Hot or Cold Temperature

 Do Not Resterilize

 Do Not Reuse

 Do Not Use if Package is Damaged

 Keep Dry

 Manufacturer

 MR Safe

 Serial Number

[STERILE R] Sterilized using Irradiation

 Use By

 Configured For Hiatal Hernia Repair

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© 2007-2011, 2014-2015, 2018-2020 W. L. Gore & Associates, Inc.
GORE®, BIO-A®, and designs are trademarks of W. L. Gore & Associates.
MINIPAX® is a trademark of Multisorb Technologies, Inc.
W. L. GORE & ASSOCIATES B.V.