INSTRUCTIONS FOR USE FOR:

en
English **hu**
Magyar

cs
Čeština **it**
Italiano

da
Dansk **no**
Norsk

nl
Nederlands **pl**
Polski

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fr
Français **es**
Español

de
Deutsch **sv**
Svenska

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
INSTRUCTIONS FOR USE FOR GORE® SYNECOR INTRAPERITONEAL BIOMATERIAL

Carefully read all instructions prior to use. Observe all instructions, warnings, and precautions noted throughout. Failure to do so may result in complications.

PRODUCT DESCRIPTION

The GORE® SYNECOR Intraperitoneal Biomaterial is a composite mesh comprised of a non-absorbable macroporous knit constructed of monofilament polytetrafluoroethylene (PTFE) fibers. The visceral surface is a bioabsorbable synthetic nonporous film comprised of poly (glycolide:trimethylene carbonate) copolymer (PGA:TMC). The parietal surface is a bioabsorbable synthetic porous fibrous structure comprised of PGA:TMC copolymer. Degraded via a combination of hydrolytic and enzymatic pathways, the PGA:TMC copolymer has been found to be both biocompatible and non-immunogenic. In vivo studies with this copolymer indicate the bioabsorption process should be complete by six to seven months.¹ GORE® SYNECOR Intraperitoneal Biomaterial is designed for intraperitoneal placement.

INDICATIONS FOR USE

The GORE® SYNECOR Intraperitoneal Biomaterial device is intended for use in the repair of ventral/incisional hernias that may require the addition of a non-absorbable reinforcing or bridging material.

CONTRAINDICATIONS

The GORE® SYNECOR Intraperitoneal Biomaterial is contraindicated for use in reconstruction of cardiovascular defects.

WARNINGS AND PRECAUTIONS

Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient (see ADVERSE EVENTS).

WARNINGS

GENERAL

- Do not use if the package is opened or damaged, or it is suspected that the sterility of the device has been compromised, as infection and related serious potential patient harms could occur.
- Do not use past the "use by" (expiration) date, as this may affect product performance and sterility, resulting in serious potential harms.
- Reuse may result in infection and related serious potential patient harms. See HOW SUPPLIED.
- Use of this product in applications other than those indicated, including but not limited to those listed below, has the potential for serious complications:
 - The safety and effectiveness of this device has not been established for use in pelvic floor soft tissue deficiencies. Clinical data has shown that use of any surgical mesh in pelvic organ prolapse and stress urinary incontinence procedures could result in complications, including, but not limited to, mesh exposure, mesh erosion, pelvic pain, fistula, dyspareunia, infection, vaginal bleeding or discharge, vaginal dysfunction or recurrence, which may require additional surgery. Further, a positive risk benefit profile for surgical mesh used in transvaginal pelvic organ prolapse repair is not well established in published literature.
 - The safety and effectiveness of this device has not been established for use in patients who have a surgical site infection that cannot be treated successfully prior to device placement (e.g. use of surgical mesh in an infected surgical site). Use of any surgical mesh in the presence of infection could result in potentially serious patient harms and additional intervention, including surgery.

BEFORE IMPLANT PROCEDURE

- As with any mesh, use of GORE® SYNECOR Intraperitoneal Biomaterial in the presence of contamination may result in fever, infection, irritation or inflammation, pain, and wound dehiscence and may require removal of the mesh if infection occurs.

PROCEDURE

- Strict aseptic techniques should be followed. Failure to do so may result in infection and related serious potential patient harms.
- Failure to ensure adequate device overlap of the defect may result in recurrence of the original, or development of an adjacent tissue defect which may lead to bowel obstruction, pain, and secondary procedure.
- As GORE® SYNECOR Intraperitoneal Biomaterial is designed for intraperitoneal placement, the shiny film surface should not be placed adjacent to fascial or subcutaneous tissue because it may result in minimal tissue attachment. Persistent seroma may result and may require additional intervention including surgery to resolve.
- Correct surface orientation is important for GORE® SYNECOR Intraperitoneal Biomaterial. Incorrect parietal (non-shiny) surface placement with unnecessary contact with internal organs may result in adhesions, erosion, extrusion, fistula, infection, irritation or inflammation, pain, and additional intervention including surgery. Clinical data on abdominal surgery, with or without surgical mesh, has shown that adhesion formation may occur and may result in complications including pain, bowel obstruction, and additional intervention including surgery.

- Use of absorbable sutures to secure the mesh in place may result in inadequate fixation. Ensure staple size and staple or tack spacing provides adequate fixation of device to tissue. Insufficient fixation of the device may result in mesh migration, erosion or extrusion, or hernia recurrence, which may lead to bowel obstruction, fistula, pain, and additional intervention including surgery.

AFTER IMPLANT PROCEDURE

- When post-operative infection is suspected, debridement of involved tissues should be considered. If an infection develops, it should be treated aggressively. An unresolved infection may require removal of the material.
- If using this device as a permanent implant and exposure occurs, treat to avoid contamination, or device removal may be necessary. Exposure or protrusion of the device could result in infection, pain, and additional intervention including surgery.

PRECAUTIONS

The following caution should be considered to ensure the safe and effective use of the device and protect the safety of the patient.

- Use of GORE® SYNECOR Intraperitoneal Biomaterial in infants, children or pregnancy where future growth will occur has not been evaluated.

ADVERSE EVENTS

POTENTIAL CLINICAL AND DEVICE ADVERSE EVENTS

Device-related:

Possible adverse reactions with the use of any tissue deficiency prosthesis may include, but are not limited to adhesions and related harms, bowel obstruction, contamination, defect recurrence and related harms, erosion or extrusion and related harms, exposure or protrusion and related harms, fever, fistula, infection, irritation or inflammation, mesh migration, mesh shrinkage, pain, paresthesia, seroma or hematoma and related harms, tissue ischemia, wound dehiscence, and additional intervention including surgery.

Procedure-related:

As with any surgical procedure, there are always risks of complications for surgical repair of hernias, with or without mesh, these may include but are not limited to, adhesions and related harms, bleeding, bowel obstruction, erosion or extrusion and related harms, exposure or protrusion and related harms, fever, fistula, hernia recurrence, ileus, increased procedure time and related harms, irritation or inflammation, infection, pain, paresthesia, perforation, revision/resurgery, seroma or hematoma and related harms, wound complications and wound dehiscence.

Refer to IFU Section WARNING AND PRECAUTIONS for any additional information regarding adverse events and any steps that should be taken to avoid them, as well as information about other warnings and precautions.

CLINICAL AND DEVICE ADVERSE EVENT REPORTING

Any adverse event involving the GORE® SYNECOR Intraperitoneal Biomaterial should be reported to the manufacturer and the country specific regulatory authorities immediately. To report an event to W. L. Gore & Associates, email: medcomplaints@wlgore.com or contact:

US:

Phone: +1.800.528.1866 or +1.928.864.4922,

Fax: +1.928.864.4364

EMEA:

Phone: +49 89 4612 3440,

Fax: +49 89 4612 43440

To report an adverse event to W. L. Gore & Associates in all other markets, please use the email above (medcomplaints@wlgore.com) or call the US phone number (+1.800.528.1866 or +1.928.864.4922).

HOW SUPPLIED

CONTENTS

The GORE® SYNECOR Intraperitoneal Biomaterial is packaged in a single sterile barrier foil pouch with an internal TYVEK® desiccant pocket. The pouch is then placed in an envelope style paperboard carton.

STERILITY

- The GORE® SYNECOR Intraperitoneal Biomaterial is provided STERILE and non-pyrogenic, and is sterilized by ethylene oxide. Do not resterilize.
- Do not use after the "use by" (expiration) date printed on the label. Provided that the integrity of the package is not compromised in any way, the package will serve as an effective barrier until the "use by" (expiration) date printed on the label.
- The GORE® SYNECOR Intraperitoneal Biomaterial is designed for single use only; do not reuse device. Gore does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination which may lead to patient harms. (See WARNINGS)

STORAGE AND HANDLING

- Store in a dry place. Avoid exposing the package or device to extreme hot or cold temperatures.
- Do not use the GORE® SYNECOR Intraperitoneal if the foil pouch is compromised or if the device is damaged.
- To open the package, peel the seal apart from the upper portion of the foil pouch.
- Use clean sterile gloves oratraumatic sterile instruments when handling the GORE® SYNECOR Intraperitoneal Biomaterial.
- Handle and dispose of the device and packaging, taking into account any infectious or microbial hazard risks, with the necessary precautionary measures in accordance with acceptable medical practice and with applicable local, state, and federal laws and regulations.
- See WARNINGS AND PRECAUTIONS for additional considerations specific to storage and handling.

For additional Storage and Handling information, see HOW SUPPLIED.

DIRECTIONS FOR USE

PHYSICIAN TRAINING REQUIREMENTS

The GORE® SYNECOR Intrapерitoneal Biomaterial should only be used by physicians experienced in the repair of hernias.

DEVICE PREPARATION

- Check information on the product labeling for correct configuration, size, and expiration date.
- Open the carton at the flap on the top of the package and remove the pouch.
- Inspect the package and do not use device if the package is damaged or if it is suspected sterility has been compromised.
- Use sterile technique to peel the seal apart from the upper portion of the pouch and remove the device.
- Ensure the size of the mesh is adequate for the intended repair. Use sharp surgical instruments to trim GORE® SYNECOR Intrapерitoneal Biomaterial to ensure appropriate overlap of the defect on all sides. If GORE® SYNECOR Intrapерitoneal Biomaterial is cut too small, excessive tension may be placed on the suture line, which may result in recurrence of the original, or development of an adjacent, tissue defect. SEE WARNINGS.
- GORE® SYNECOR Intrapерitoneal Biomaterial does not require pre-wetting. If desired, GORE® SYNECOR Intrapерitoneal Biomaterial may be wetted with sterile saline or similar standard irrigation solution prior to placement through a trocar to facilitate introduction.

See WARNINGS AND PRECAUTIONS for additional considerations specific to the procedure.

IMPLANT PROCEDURE

DEVICE INSERTION

- Insert the device through an open incision, trocar incision, or introduce through a trocar. Ease of trocar introduction may vary depending on rolled device size and trocar availability. Excessive force should be avoided during trocar introduction to minimize potential for device damage.

SURFACE ORIENTATION

- Ensure correct surface orientation during intraperitoneal placement. One surface of the product has a shiny nonporous film. The shiny film surface should be placed adjacent to those tissues where minimal tissue attachment is desired (i.e., serosal surfaces). The parietal (non-shiny) surface should be placed adjacent to those tissues where ingrowth is desired. SEE WARNINGS.

SUTURE FIXATION

- Nonabsorbable sutures, such as GORE-TEX® Suture, of appropriate size with a taper or piercing point needle are recommended to secure the mesh. The use of absorbable sutures may lead to inadequate anchoring of GORE® SYNECOR Intrapерitoneal Biomaterial to the host tissue and necessitate reoperation. SEE WARNINGS.
- For best results, use monofilament sutures. Suture size should be determined by surgeon preference and the nature of the reconstruction.
- When suturing GORE® SYNECOR Intrapерitoneal Biomaterial to the host tissue, the mesh should be sutured a minimum distance of 1 cm (0.4 in) from the edge of the device.
- Interrupted sutures can provide additional security against recurrence due to suture failure.

ALTERNATIVE FIXATION

- Permanent staples or helical tacks (also known as helical coils) can be used as an alternative to sutures.
- Use of absorbable tack and staple fixation devices with GORE® SYNECOR Intrapерitoneal Biomaterial have been evaluated acutely in animal models and shown to be compatible.
- Appropriate force should be used when positioning tacking devices to avoid damage to film surface.
- Staple size and staple or tack spacing should be determined by surgeon preference to provide for adequate tissue fixation and to prevent reherniation. SEE WARNINGS.
- Fixation methods other than those described above have not been evaluated for use with GORE® SYNECOR Intrapерitoneal Biomaterial.

All fixation devices should be used in accordance with the manufacturer's instructions for use. The physician should reference the manufacturer's instructions for use and any supporting clinical literature regarding any potential warnings, precautions, and adverse reactions related to fixation devices.

See WARNINGS AND PRECAUTIONS for additional considerations specific to the procedure.

MR SAFETY INFORMATION

GORE® SYNECOR Intrapерitoneal Biomaterial is MR Safe.

REFERENCES

1. Katz AR, Mukherjee DP, Kaganov AL, Gordon S. A new synthetic monofilament absorbable suture made from polytrimethylene carbonate. *Surgery, Gynecology & Obstetrics* 1985;161(3):213-222.

DEFINITIONS

[EC] REP Authorised Representative in the European Community

[REF] Catalogue Number

 Caution

RX Only CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

 Consult Instructions for Use

 Date of Manufacture

 Do Not Expose to Extreme Hot or Cold Temperature

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© 2016, 2017, 2020 W. L. Gore & Associates, Inc.
GORE®, GORE-TEX®, SYNECOR, and designs are trademarks of W. L. Gore & Associates.
W. L. GORE & ASSOCIATES B.V.