



Breathe easy

Respironics AF531 mask with Aerogen NIVO nebulizer system

With the Respironics AF531 oro-nasal mask and Aerogen NIVO nebulizer system (Pro-X controller, NIVO nebulizer, and NIVO EE or SE elbow) clinicians can provide nebulizer treatments without disconnecting patient circuits and struggling with in-line nebulizer add-ons. The nebulizing elbows are designed specifically for noninvasive ventilation (NIV) and can remain attached when not in use or be easily removed between treatments.

Key advantages

- Supports continuous NIV and medication delivery
- Delivers medication in upright and reclining positions
- Uses Aerogen OnQ™ vibrating mesh technology and Pro-X controller

PHILIPS
sense and simplicity

Respironics AF531 with Aerogen NIVO nebulizer system

The AF531 can be used with the nebulizing NIVO leak 1 entrainment (EE) or NIVO standard elbow (SE) and a Aerogen NIVO nebulizer.



Respironics AF531 mask with nebulizer elbow

Description	Quantity	Part number
AF531 with four-point headgear and NIVO leak 1 entrainment elbow (EE), small	10	1072633
AF531 with four-point headgear and NIVO leak 1 entrainment elbow (EE), medium	10	1072634
AF531 with four-point headgear and NIVO leak 1 entrainment elbow (EE), large	10	1072635
AF531 with CapStrap headgear and NIVO leak 1 entrainment elbow (EE), small	10	1072636
AF531 with CapStrap headgear and NIVO leak 1 entrainment elbow (EE), medium	10	1072637
AF531 with CapStrap headgear and NIVO leak 1 entrainment elbow (EE), large	10	1072638
AF531 procedure pack with circuit and NIVO leak 1 entrainment elbow (EE), small	5	1072642
AF531 procedure pack with circuit and NIVO leak 1 entrainment elbow (EE), medium	5	1072643
AF531 procedure pack with circuit and NIVO leak 1 entrainment elbow (EE), large	5	1072644



NIVO leak 1
entrainment elbow (EE)



NIVO standard elbow (SE)



NIVO nebulizers

Description	Quantity	Part number
NIVO nebulizer	5	1076302
NIVO nebulizer and NIVO standard elbow (SE)	5	1076303
NIVO nebulizer and NIVO leak 1 entrainment elbow (EE)	5	1076304

NIVO replacement elbows

Description	Quantity	Part number
NIVO leak 1 entrainment elbow (EE)	5	1072650
NIVO leak 1 entrainment elbow (EE)	1	1072649
NIVO standard elbow (SE)	5	1076754
NIVO standard elbow (SE)	1	1076755

Nebulizer systems and Aerogen controller kits

NIVO Pro-X nebulizer systems*

Description	Part number
NIVO Pro-X nebulizer system - English language EU	1061733
NIVO Pro-X nebulizer system - Northern EU	1061734
NIVO Pro-X nebulizer system - Southern EU	1061735
NIVO Pro-X nebulizer system - Scandanavian EU	1082689
NIVO Pro-X nebulizer system - UK and Ireland	1072619
NIVO Pro-X nebulizer system - Australia and New Zealand	1072620
NIVO Pro-X nebulizer system - English language international	1072633
NIVO Pro-X nebulizer system - International	1072634
NIVO Pro-X nebulizer system - Japan	1072635

*All NIVO Pro-X nebulizer systems include one Pro-X controller; one power adapter; one mounting bracket, one mounting clamp, one cable, one OIS, one NIVO leak 1 entrainment elbow (EE), one NIVO standard elbow (SE), and 2 NIVO aerosol generators.



NIVO Pro-X controller kits*

Description	Part number
NIVO Pro-X controller kit - English language EU	1079421
NIVO Pro-X controller kit - Northern EU	1082693
NIVO Pro-X controller kit - Southern EU	1082694
NIVO Pro-X controller kit - Scandanavian EU	1082695
NIVO Pro-X controller kit - UK and Ireland	1079422
NIVO Pro-X controller kit - Australia and New Zealand	1079423
NIVO Pro-X controller kit - English language international	1079424
NIVO Pro-X controller kit - International	1082696
NIVO Pro-X controller kit - Japan	1082698

*All NIVO Pro-X controller kits include one Pro-X controller; one power adapter; one mounting bracket, one mounting clamp, one cable, and one OIS.

Respironics has partnered with Aerogen to develop a nebulizing elbow that interfaces with the AF531 full face mask and the Aerogen NIVO nebulizer and Pro-X controller.

Please visit www.philips.com/nivo



© 2011 Koninklijke Philips Electronics N.V. All rights are reserved.
Aerogen, Aeroneb, and NIVO are trademarks of Aerogen, Inc. CapStrap and Respironics is a trademark of Respironics, Inc. and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 71241 * APR 2011