



# Optimizer™ Polyp Trap

## Multiple Chamber Polyp Trap

The versatility of the Optimizer™ Polyp Trap is highlighted by our independent, selectable specimen chambers, four direct bypass channels and a convenient central vacuum lumen. A direct connection to wall suction makes for fast and easy fluid management.

Each trap is individually packaged for single-patient use to reduce risk of **cross-contamination**.



## Features and Benefits

### Versatile

- Four independent, selectable trap chambers
- Multiple direct bypass channels
- Visual markers for each chamber

### Fast

- Convenient central vacuum lumen

### Individually Packaged

- Single patient use reduces risk of cross-contamination



Direct connection to wall suction for easy fluid management



Sequenced trap chambers and direct bypass chambers

## Ordering Information

Description	Qty	Catalog Number
Optimizer™ Polyp Trap Multi Chamber Polyp Trap (Non-Sterile, Disposable)	20/Box	CAP132

CONMED Corporation  
525 French Road, Utica, NY 13502

Local 727-392-6464  
Toll Free 1-800-237-0169

CONMED.com  
customer\_service@CONMED.com