

ScopeSaver™ Latex-Free Bite Block



It's All In The Strap

93% of G.I. Nurses Chose the ConMed Latex-Free Strap for Patient Comfort When Compared to the leading Competitor*

Applied Ingenuity™

 **CONMED™**
ENDOSCOPIC TECHNOLOGIES

ScopeSaver™ Latex-Free Bite Block



DESIGNED FOR PROTECTION AND COMFORT

- Our straps are designed to reduce entanglement during application and removal of bite block.
- Smooth Strap for patient comfort
- Protect your staff and patients from Latex Exposure
- Prevent costly scope repairs

SUPERIOR BITE BLOCK DESIGN

- Large, flexible side port openings allow insertion of suction tube and enables finger-assisted scope passage
- Rear-guard retaining lip maximizes Bite-Block retention and scope protection
- The large Scope Saver Accommodates up to 60F American Dilator

*Data on file

Scope Saver® Disposable Bite Blocks

Latex **FREE**

Non-Sterile

Catalog Number	Description	Lumen Size	Color	Unit
000429*	Large Scope Saver® Disposable Latex Free Bite Block	20 x 27mm	Green	50/Box
100429*	Small Scope Saver® Disposable Latex Free Bite Block	14 x 20mm	Purple	25/Box

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



ConMed Endoscopic Technologies
129 Concord Road
Billerica, MA 01821

Customer Service: 1-800-448-6506
FAX: 1-800-438-3051
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
email: info@mail.conmed.com
www.conmed.com

© 2010, CONMED Corporation, 10/2010, Control #12617 Rev C

MBR 9131 Rev 2