

# RELIANCE

## Bulldog Clamps



INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*

—

Reliance Bulldog Clamps


—

Reliable, Surgeon-Controlled Vessel Clamps

The Scanlan® Reliance Bulldog Clamp has been designed to provide secure, robust clamping pressure to assist in the temporary occlusion of vessels. Manufactured with the highest quality stainless steel, its 1 x 2 (25 mm) DeBakey jaws ensure a firm yet atraumatic clamping pressure. A unique suture port allows optional suture tagging to enhance surgical location and recovery.

Features Include...

- 1 x 2 (25 mm) DeBakey jaws
- 350 grams closing pressure\*
- Suture port
- Highest quality stainless steel
- Reusable / durable
- Low profile design
- Radiopaque



Reliance Bulldog Applier / Introducer

- Custom jaw allows for secure deployment of bulldog into the surgical field
- Designed for use with 12 mm cannula



\* Approximate clamping pressure

—

Reliance Bulldog Clamps

—

Reliable, Surgeon-Controlled Vessel Clamps

Bulldog Clamps



**3795-50** (1 unit)\*\*  
25 mm straight DeBakey jaws,  
2" (5 cm) 350 grams closing pressure\*



**3795-51** (1 unit)\*\*  
25 mm curved DeBakey jaws,  
2" (5 cm) 350 grams closing pressure\*



**3795-54** (1 unit)\*\*  
7 mm straight DeBakey jaws,  
1 1/4" (3.2 cm) 350 grams closing pressure\*  
*Ideal for selective regional clamping*

**9999-50** (4 units plus tray)\*\*  
Starter Pack - Includes 4 Bulldog Clamps  
(3795-50) and Metal Sterilization Tray (2133-03)

\*\* Firm Sale

Applier / Introducer  
9009-395

Working length: 13 3/4" (35 cm)  
Total length: 17 3/4" (44 cm)  
Shaft diameter: 7 mm

# — Reliance Bulldog Clamps —

Reliable, Surgeon-Controlled Vessel Clamps

## Metal Sterilization Tray



### 2133-03

Micro Sterilization Tray with mat and hardware  
(Reliance Bulldog Clamps not included)

**For more information or to place an order,**  
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458  
International: 651-298-0997  
FAX: 651-298-0018  
[www.scanlaninternational.com](http://www.scanlaninternational.com)



INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*  
ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.  
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • [www.scanlaninternational.com](http://www.scanlaninternational.com)

MK-B-RBDC Rev. E  
© 2012 Scanlan International, Inc.