




VENCLOSE™

# RF Ablation System



## Streamlining the Treatment of Venous Reflux Disease

The next-generation endovenous radiofrequency ablation solution

### VENCLOSE™ RF Ablation Catheter

#### Designed for Versatility

- » Two heating lengths (10 cm and 2.5 cm) in one catheter reduces total in-vein time by >30%\* and simplifies inventory management
- » Smaller, 6 Fr profile minimizes invasiveness


#### Optimized for Efficiency

- » Markings and curved tip ease navigation and positioning
- » Smaller profile offers greater flexibility

### VENCLOSE™ RF Generator

#### Simplified Setup and Operation



- » Touchscreen display shows real-time procedure data to inform treatment decisions
- » Simple catheter connection port allows quick and easy 'blind' catheter plug-in
- » Nearly instant power-up
- » Compact design
- » Audible tones for thermal delivery allows you to focus on the patient, not the display



\*When compared to the ClosureFast™ Procedure. ClosureFast is a trademark of a Medtronic company.

## Efficient Thermal Delivery

Built on the well-documented history of coil-based endovenous RF ablation.



## Ordering Information

Description	Order Number
VENCLOSE™ RF Generator	VC-RFG-1
VENCLOSE™ RF Ablation Catheter (60 cm)	VC-10A25-6F-60
VENCLOSE™ RF Ablation Catheter (100 cm)	VC-10A25-6F-100
Optional foot pedal	VC-FP-1
Region-specific power cord (where X is region-dependent)	VC-PC-X



VENCLOSE, Inc.  
2570 N. First Street, 2nd Floor #221  
San Jose, CA 95131, USA

Phone: +1-844-834-6292  
info@venclose.com  
**venclose.com**

<sup>1</sup>. Data on file.  
Copyright © VENCLOSE, Inc. All rights reserved.  
VENCLOSE is a trademark of VENCLOSE in the United States and other countries.