



Clarius HD3 Vet Scanners

Clear Ultrasound Imaging
You Can Trust



- ✓ High-Definition Ultrasound
- ✓ Steadfast Wireless Connection
- ✓ 3-Year Standard Warranty
- ✓ Unlimited Users
- ✓ 1-Hour Remote In-Service Training
- ✓ M-Mode, Power Doppler, Color Doppler
- ✓ Easy-To-Use App for iOS and Android
- ✓ Specialized Scanners for Optimal Imaging

Superior Image Quality

- High-quality images that are comparable to cart-based systems
- Imaging settings are optimized based on scanning application
- Optimal results from specialized scanners

Reliable Portability

- Wireless freedom: no cables to trip over or break, and less risk of infections during procedures
- Stable connectivity from our proprietary multi-channel Wi-Fi architecture

Easy to Use

- Point-and-Shoot Ultrasound™ technology removes the need to adjust knobs and switches
- Artificial intelligence adjusts parameters in real-time
- Compact and waterproof device that feels like standard scanners

Model	Frequency	Max Depth	Specialties
C3 HD3 Vet	2-6 MHz	40 cm	Large Animals (Large Canine, Goat, Sheep, Cows, Equine)
C7 HD3 Vet	3-10 MHz	18 cm	Small-Medium Animals (Cats, Dogs)
L7 HD3 Vet	4-13 MHz	11 cm	Small-Medium Animals, Superficial Studies (Equine MSK)