





Designed for exceptional  kink resistant performance and non-ovaling effect.

# DuraSheath

Kink resistant introducer sheath

 Rethink the basics


# Supporting complex vascular access

## DuraSheath

### Kink resistant introducer sheath

From easy insertion to smooth navigation to superb visibility: DuraSheath is designed to combine optimal performance with excellent pushability and trackability.

While maintaining the inner diameter, DuraSheath allows for smooth device passage with an exceptionally small kink radius and non-ovaling effect. The detachable hemostatic valve provides flexibility and versatility, allowing you the ability to focus on the intervention.



#### **New optional shaft material**

optimal backup support transitioning to softer distal zone for safe navigation.



#### **Dilator with snap-fit**


secures the dilator for safe insertion.



#### **Detachable hemostatic valve**

provides you options.





**Hydrophobic coating**  
allowing for effortless & durable gliding  
properties, no matter if wet or dry.

## 🎯 Rethink the basics

Passionate about bringing innovative ideas and medical device solutions to life, we are a team building on nearly 20 years of experience in the design and manufacturing of medical device solutions for partners around the world.

Even small changes in product design and quality can make a clear difference – making procedures easier and simpler with better outcomes.

Continually striving for solutions to the most basic problems, we don't go back to basics, we rethink the basics.



**Atraumatic tip**  
designed for smooth and easy insertion.




**Platinum iridium  
marker band**  
for excellent radiopacity.



Shaft versions  
color coded.

**new** Optional tip designs


**MP** Multipurpose




**LIMA** Left Internal Mammary Artery




**RDC** Renal Double Curve



**HS** Hockey Stick



**Outer diameter**  
smooth surface for  
low friction navigation.




**Inner diameter**  
smooth device delivery  
and movement.

# DuraSheath

Technical details / Ordering information

Providing access to target vessel not limited to

- peripheral
- renal
- coronary
- carotid



Access all areas

## Introducer Sheath

Material – outer	Stainless steel reinforced polymer
Material – inner	PTFE liner
Flexibility	Regular & Firm
Sheath coating	Hydrophobic
Sheath tip	Smooth taper with radiopaque platinum iridium marker band
French size	4, 5, 6, 7, 8 Fr
Length	15, 45, 60, 90, 130 cm
Sidearm/hemostatic valve	Detachable, with color coded 3-way-stopcock

Tip design									
		Regular	Firm	Firm	Firm	Firm	Firm		
Diam. (Fr)	Length (cm)	Catalogue Number						Inner Diameter Min. (mm/inch)	Units/Case
4 Fr	15 cm	FG-04015						1.60 mm / 0.063"	5
4 Fr	45 cm	FG-04045						1.60 mm / 0.063"	5
4 Fr	60 cm	FG-04060						1.60 mm / 0.063"	5
4 Fr	90 cm	FG-04090						1.60 mm / 0.063"	5
4 Fr	130 cm		FG-04130-F					1.60 mm / 0.063"	5
5 Fr	15 cm	FG-05015						1.92 mm / 0.075"	5
5 Fr	45 cm	FG-05045	FG-05045-F	FG-05045-F-MP	FG-05045-F-LIMA	FG-05045-F-RDC	FG-05045-F-HS	1.92 mm / 0.075"	5
5 Fr	60 cm	FG-05060	FG-05060-F					1.92 mm / 0.075"	5
5 Fr	90 cm	FG-05090	FG-05090-F	FG-05090-F-MP				1.92 mm / 0.075"	5
6 Fr	15 cm	FG-06015	FG-06015-F					2.25 mm / 0.088"	5
6 Fr	45 cm	FG-06045	FG-06045-F	FG-06045-F-MP	FG-06045-F-LIMA	FG-06045-F-RDC	FG-06045-F-HS	2.25 mm / 0.088"	5
6 Fr	60 cm	FG-06060	FG-06060-F					2.25 mm / 0.088"	5
6 Fr	90 cm	FG-06090	FG-06090-F	FG-06090-F-MP				2.25 mm / 0.088"	5
6 Fr	130 cm		FG-06130-F					2.25 mm / 0.088"	5
7 Fr	15 cm	FG-07015						2.55 mm / 0.100"	5
7 Fr	45 cm	FG-07045	FG-07045-F	FG-07045-F-MP	FG-07045-F-LIMA	FG-07045-F-RDC	FG-07045-F-HS	2.55 mm / 0.100"	5
7 Fr	60 cm	FG-07060	FG-07060-F					2.55 mm / 0.100"	5
7 Fr	90 cm	FG-07090	FG-07090-F	FG-07090-F-MP				2.55 mm / 0.100"	5
8 Fr	15 cm	FG-08015						2.88 mm / 0.113"	5
8 Fr	45 cm	FG-08045	FG-08045-F					2.88 mm / 0.113"	5
8 Fr	60 cm	FG-08060	FG-08060-F					2.88 mm / 0.113"	5
8 Fr	90 cm	FG-08090	FG-08090-F					2.88 mm / 0.113"	5

15 cm versions contain two dilators.

## Dilator Guidewire Compatibility

Diameter	4 Fr / 5 Fr / 6 Fr		7 Fr / 8 Fr	
Length	15 cm	45 cm / 60 cm / 90 cm / 130 cm	15 cm	45 cm / 60 cm / 90 cm
Dilator Guidewire Comp. max.	0.035" & 0.018"		0.035"	0.035" & 0.018" / 0.038"

**CMI**  A product of  
Contract Medical International GmbH  
Lauensteiner Str. 37  
01277 Dresden • Germany  
Tel. +49 351 213 88 88  
[HeraeusMedicalComponents.com](http://HeraeusMedicalComponents.com)

Contract Medical International  
joins Heraeus Medical Components

**CMI**  | **Heraeus**