

PHILIPS

MR Patient Care

Philips MR Patient Monitoring solutions

Premium MR Patient Monitor for high-acuity care Expression MR400



- 15" LED touchscreen
- 5,000 Gauss
- 4 W/kg SAR & 7.2 μT B1_{rms}
- 1.5T & 3.0T MR compliance
- Gauss alarm (with AA option)
- Parameters:
 - Basic: ECG*, SpO₂*, NIBP, CO₂/Respiration
 - Optional: IBP (x2), body temperature, anesthetic agents
- Advanced clinical alarms
- MR gating feature: wired/wireless

* wireless ECG & SpO₂

Premium MR Patient Monitor for low-acuity care Expression MR200



- 12" LCD
- 1,500 Gauss
- 4 W/kg SAR
- 1.5T & 3.0T MR compliance
- Gauss alarm (standard)
- Parameters:
 - Basic: ECG*, SpO₂*, NIBP
 - Optional: CO₂/respiration
- MR gating feature: wired/wireless

* wireless ECG & SpO₂



Information Portal MR Patient Care Portal 5000

- Medical grade projected capacity touch panel 18.3" LED display
- Remote monitoring with Expression MR200 and Expression MR400
- Modular design for cost effective and possible future upgrades
- HIS/EMR connectivity with enhanced product and IT security
- Patient management
- Password protected administrative profile features
- Advanced visual perception principles and decision-making tools
- Optional: barcode scanner, printer, wireless keyboard and mouse

More than 30 years of excellence in MR patient monitoring

Key facts at a glance

Patient types

- Adults
- Infants
- Pediatric
- Neonates



Parameters & alarms

- Bedside-level
- Advanced clinical alarms



Make the most of the MR

4 W/kg SAR | 5,000 Gauss*
7.2 μT B1_{rms}* | 1.5T & 3.0T

Connect

- with convenience
- Enhanced product and IT security
- Wireless communication
- HL7/RS232 connectivity
- HIS/EMR



* Based on Expression MR400 specifications

Your Philips contact:

