


DIGITAL VIDEO MONITOR

1280 X 800 PIXELS


Your Trusted ICU Partner

LONG STANDBY TIME PORTABLE MONITOR



VATHiN

Version 20005



VATHIN Vathin Medical

SPECIFICATION / DIGITAL VIDEO MONITOR /

ELECTRICAL POWER

Medical Adapter	100V~240VAC 50/60Hz 1.4~0.7A
Power Requirement	15 V, 4.0 A DC 60W input
Battery Type	10.8 V 10.05Ah Li-ion
Battery Operation	More than 8 hours

DISPLAYPLAY

Max. Resolution	1280 x 800 Pixels
Orientation	Gyro
Display Type	12.8" Colour TFT LCD
Image Adjustment	Contrast and Sharpness , Brightness

FUNCTION

Snapshot	1089x800p /Png
Video Record	720x576p /MP4
AWB	Yes
Playback	Picture & Video

CONNECTIVITY

USB Connection	2 x USB 3.0 Type A
Digital Video out	HDMI (1280 x 800p, 60 fps)
WiFi	Supports IEEE 802.11ac/a/b/g/n
LAN	Ethernet RJ45 connector, 10/100/1000Mb

DIMENSIONS

Width x Height x Thickness: 308 mm x230mmx43m
Weight : 2.3 kg

Version 20005