





Reinforce with confidence



Pulmonary Resection



Gastric Resection



Sleeve Gastrectomy

PERFORMANCE through innovation


GORE® SEAMGUARD® Bioabsorbable Staple Line Reinforcement consists of a unique, patented synthetic tissue scaffold that provides consistent, predictable performance, every time.

The product is engineered to reduce the incidence of perioperative leaks and bleeding in a variety of open and minimally invasive surgeries^{1,2}. Then, after six months, GORE® SEAMGUARD® Bioabsorbable Staple Line Reinforcement is completely absorbed, eliminating the risk of a prolonged inflammatory response.

- 100% bioabsorbable, synthetic material
- Proven material
- Predictable performance
- Strong, yet thin



Intestinal Anastomosis



Left Atrial Appendage Ligation



SEAMGUARD®
BIOABSORBABLE

STAPLE LINE REINFORCEMENT

1. Franklin. Safety and efficacy of the use of bioabsorbable "Seamguard" in colorectal surgery at the Texas Endosurgery Institute. *Surgical Laparoscopy, Endoscopy & Percutaneous Techniques* 2005;15(1):9-13.

2. Consten EC, Gagner M, Pomp A, Inabnet WB. Decreased bleeding after laparoscopic sleeve gastrectomy with or without duodenal switch for morbid obesity using a stapled buttressed absorbable polymer membrane. *Obesity Surgery* 2004;14:1360-13 66.

Configured For Endoscopic Staplers

	STAPLING DEVICE	CARTRIDGE COLOR	CATALOGUE NUMBER
Ethicon Endo-Surgery	ENDOPATH® ETS 45 Endoscopic Transecting	White, Blue, Green	1BSGETS45 / 12BSGETS45*
	ENDOPATH® ETS Flex 45	White, Blue, Green	1BSGETS45 / 12BSGETS45*
	ENDOPATH® EZ 45	Blue, Green	1BSGEZ45 / 12BSGEZ45*
	ENDOPATH® ETS Flex 60	White, Blue, Green	1BSGETS60 / 12BSGETS60*
	Echelon 45 ENDOPATH® Stapler	White, Blue, Gold, Green	1BSGEC45 / 12BSGEC45*
	Echelon 60 ENDOPATH® Stapler	White, Blue, Gold, Green	1BSGEC60 / 12BSGEC60*
	Echelon FLEX 60 ENDOPATH®	White, Blue, Gold, Green	1BSGEC60A / 12BSGEC60A*
Covidien Autosuture	Endo GIA Universal Straight 45	White, Blue	1BSGUSS45B / 12BSGUSS45B*
	Endo GIA Universal Reticulator 45	White, Blue	1BSGUSS45B / 12BSGUSS45B*
	Endo GIA Universal Straight 45	Green	1BSGUSS45G / 12BSGUSS45G*
	Endo GIA Universal Reticulator 45	Green	1BSGUSS45G / 12BSGUSS45G*
	Endo GIA Universal Straight 60	White, Blue	1BSGUSS60B / 12BSGUSS60B*
	Endo GIA Universal Reticulator 60	White, Blue	1BSGUSS60B / 12BSGUSS60B*
	Endo GIA Universal Straight 60	Green	1BSGUSS60G / 12BSGUSS60G*
	Endo GIA Universal Reticulator 60	Green	1BSGUSS60G / 12BSGUSS60G*

* Twelve packages per box

Configured For Open Staplers

	STAPLING DEVICE	CATALOGUE NUMBER
Ethicon Endo-Surgery	PROXIMATE® Linear Cutter 55 (Staplers: TLC55, TCT55; Reloads: TCR55, TRT55)	1OBSGPROX55 / 12OBSGPROX55*
	PROXIMATE® Linear Cutter 75 (Staplers: TLC75, TCT75, TCD75; Reloads: TCR75, TRT75, TRD75)	1OBSGPROX75 / 12OBSGPROX75*
	PROXIMATE® Linear Cutter 100 (Staplers: TLC10, TCT10; Reloads: TCR10, TRT10)	1OBSGPROX100 / 12OBSGPROX100*
Covidien Autosuture	GIA DST 60 (Stapler: GIA6038S; Reload: GIA6038L) (Stapler: GIA6048S; Reload: GIA6048L)	1OBSGDST60 / 12OBSGDST60*
	GIA DST 80 (Stapler: GIA8038S; Reload: GIA8038L) (Stapler: GIA8048S; Reload: GIA8048L)	1OBSGDST80 / 12OBSGDST80*
	GIA DST 100 (Stapler: GIA10038S; Reload: GIA10038L) (Stapler: GIA10048S; Reload: GIA10048L)	1OBSGDST100 / 12OBSGDST100*

* Twelve packages per box

Configured For Circular Staplers

	STAPLING DEVICE	CATALOGUE NUMBER
Ethicon Endo-Surgery	ENDOPATH® 25 and Proximate 25 (Staplers: ECS25, CDH25, SDH25)	1BSGC25
	ENDOPATH® 29 and Proximate 29 (Staplers: ECS29, CDH29, SDH29)	1BSGC2829
	ENDOPATH® 33 and Proximate 33 (Staplers: ECS33, CDH33, SDH33)	1BSGC33
Covidien Autosuture	CEEA 25 (Staplers: 110214L, 111985)	1BSGC25
	CEEA 28 (Staplers: 110238L, 111987)	1BSGC2829
	CEEA 31 (Staplers: 110276L, 111989)	1BSGC31
	CEEA 34 (Stapler: 111981)	1BSGC3334
	DST Series CEEA 25 (Staplers: EEA25, EEAXL25, EEA2535, EEAXL2535, EEAORVIL25)	1BSGC25
	DST Series CEEA 28 (Staplers: EEA28, EEAXL28, EEA2835, EEAXL2835)	1BSGC2829
	DST Series CEEA 31 (Staplers: EEA31, EEAXL31)	1BSGC31
	DST Series CEEA 34 (Staplers: EEA34, EEAXL34)	1BSGC3334



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)

00800.6334.4673 (Europe)

800.437.8181 (United States)

928.779.2771 (United States)

goremedical.com

ETHICON, ENDOPATH, and PROXIMATE are trademarks of Johnson & Johnson Corporation. AUTOSUTURE and ENDO GIA are trademarks of Covidien AG or its affiliates.

Products listed may not be available in all markets. GORE®, PERFORMANCE THROUGH INNOVATION, SEAMGUARD®, and designs are trademarks of W. L. Gore & Associates.

© 2011 W. L. Gore & Associates, Inc.
AL2961-EN6 AUGUST 2011