





Product Specification


The device specification for single lobe LAA

	Code	D Distal Umbrella Diameter (mm)	D1 Proximal Cover Diameter (mm)	Size of LAmbre™ Delivery System
		LT-LAA-1622	16	22
LT-LAA-1824		18	24	
LT-LAA-2026		20	26	9F, 900mm
LT-LAA-2228		22	28	
LT-LAA-2430		24	30	
	LT-LAA-2632	26	32	10F, 900mm
	LT-LAA-2834	28	34	
	LT-LAA-3036	30	36	
	LT-LAA-3236	32	36	
	LT-LAA-3438	34	38	
	LT-LAA-3640	36	40	


The device specification for multiple lobes LAA

	Code	D Distal Umbrella Diameter (mm)	D1 Proximal Cover Diameter (mm)	Size of LAmbre™ Delivery System
		LT-LAA-1630	16	30
LT-LAA-1832		18	32	10F, 900mm
LT-LAA-2032		20	32	
LT-LAA-2234		22	34	
LT-LAA-2436		24	36	
LT-LAA-2638		26	38	


LAmbre™ LAA Delivery System with 45°x 30° Double Curves

	Code	Sheath ID (mm)	Sheath OD (mm)	Sheath Length (mm)
	SL8F45×30-900	2.7	3.5	900
	SL9F45×30-900	3.0	3.8	900
	SL10F45×30-900	3.4	4.1	900

LAmbre™ LAA Delivery System with 45° Single Curve

	Code	Sheath ID (mm)	Sheath OD (mm)	Sheath Length (mm)
	SL8F45-900	2.7	3.5	900
	SL9F45-900	3.0	3.8	900
	SL10F45-900	3.4	4.1	900

LAmbre™
LAA Closure System



Better and Safer Technology for LAA Closure

Lifetech Scientific (Shenzhen) Co., Ltd.

Add: Floor 1-5, Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China
Tel: +86-755-86026250 Fax: +86-755-86026251
Email: marketing@lifetechmed.com
Web: www.lifetechmed.com

EU authorized representative:
LifeTech Scientific (Europe) Coöperatief U.A.
Oliemolenstraat 60, 6416CB Heerlen, the Netherlands

2020-EN-1.0

CE
2797

 **Lifetech**
SCIENTIFIC



LAA Closure System

A Novel Left Atrial Appendage (LAA) Closure Device

Safe and Easy to Use

Innovative umbrella design

- Close LAA regardless of its anatomy
- Fully retrievable and repositionable

Patented anchor design for better device fixation

- Small distal hooks to engage LAA walls
- Round U-shaped Ends to trap trabeculations
- Over-sized umbrella to stent against LAA

8-10 Fr. small delivery sheath

The delivery system is available in

- double curves (45° X 30°) of distal tip
- single curve (45°) of distal tip


Global Clinical Trials Data Demonstrated LAmbré™ Safe and Easy for LAA Closure

High Implantation Success



99.6%
(n=218)
Global Study
Implantation
Success

Successful implantation regardless of LAA anatomies



Could include patients who may be rejected by other devices due to short LAA depth.

100%
(n=218)
Patient Inclusion
Regardless of
LAA Anatomies

Stable fixation

No device embolization or migration.

0.00%
(n=218)
Device
Embolization

* Implantation Success: stable placement within LAA and immediately residual flow less than 5mm.
 [1] Jai-Wun Park, Horst Sievert, et al. International Journal of Cardiology 265 (2018) 97-102.
 [2] He Huang, Yu Liu, et al. J Am Coll Cardiol Interv 2017;10:2188-94.

