



OB/NeoNatal Monitoring Supplies

PHILIPS

Changes to this edition

New product added

Multi Patient Disposable Cuffs – Multi Care

989803183351	Adult XL (navy blue with edge striping)	22
--------------	---	----

Single Care Disposable Cuffs

989803182291	Small Adult (light blue)	24
989803182301	Adult (navy blue)	24
989803182311	Adult XL (navy blue with edge striping)	24
989803182321	Large Adult (burgundy)	24

Product revised

Multi Patient Comfort Cuffs – Single Hose

M1573A	Small Adult, light blue tubing	19
M1574A	Adult, navy blue tubing	19
M1575A	Large Adult, burgundy tubing	19
M1576A	Thigh, grey tubing	19

Multi Patient Easy Care Cuffs – Single Hose

M4554B / M4554B5	Small Adult, light blue	20
M4555B / M4555B5	Adult, navy blue	20
M4556B / M4556B5	Adult Extra-Long, navy blue	20
M4557B / M4557B5	Large Adult, burgundy	20
M4558B / M4558B5	Large Adult Extra-Long, burgundy	20
M4559B / M4559B5	Thigh, grey	20

Multi Patient NIBP Cuff Kits – Single Hose

864288	Easy Care Adult Kit – 4 sizes	21
864290	Easy Care Assortment Kit – 6 sizes	21
864291	Easy Care Adult Long Kit – 6 sizes	21
M1578A	Multi-Patient Comfort Cuff Assortment Kit	21

Multi Patient Disposable Cuffs – Multi Care

989803183331	Small Adult (light blue)	22
989803183341	Adult (navy blue)	22
989803183361	Large Adult (burgundy)	22
989803183371	Thigh (grey)	22

Gentle Care Disposable Cuffs

M4574B	Small Adult, light blue	23
M4575B	Adult, navy blue	23
M4576B	Adult Extra-Long, navy blue	23
M4577B	Large Adult, burgundy	23
M4578B	Large Adult Extra-Long, burgundy	23
M4579B	Thigh, grey	23

Adult NIBP Pressure Air Hose

M1598B	Adult pressure air hose 1.5m	25
M1599B	Adult pressure air hose 3.0m	25

Philips Reusable SpO₂ Sensors

M1191B	Adult Finger Glove Sensor	26
M1191BL	Adult Finger Glove Sensor	26
M1191T	Adult Finger Glove Sensor	26
M1192A	Pediatric/Small-Adult Glove Sensor	26
M1192T	Pediatric/Small-Adult Glove Sensor	26
M1194A	Pediatric/Adult Ear Clip Sensor	27
M1196A	Pediatric/Adult Finger Clip Sensor	27
M1196S	Reusable Clip Adult SpO ₂ Sensor	27
M1196T	Pediatric/Adult Finger Clip Sensor	27
M1627A	Wristband	27

Philips Disposable SpO₂ Sensors

M1131A	Adult/Pediatric Finger Clip Disposable Sensor	28
M1133A	Neonatal/Infant/Adult Disposable Wrap Sensor (single patient use)	28
M1134A	Adhesive-Free Neonatal/Infant/Adult Disposable Wrap Sensor (single patient use)	28

Adhesive SpO₂ Sensors

M1904B	Adult Disposable Sensor	28
--------	-------------------------	----

Philips SpO₂ Adapter Cables

M1941A	SpO ₂ Extension Cable, 6.6 feet (2m)	29
--------	---	----

Single Patient Disposable NIBP Cuffs

M1820-60020	Neonatal Single Patient NIBP Cuff Kit (Sizes 1-4)	55
M1866A	Neonatal Single Patient NIBP Cuff, Size 1	55
M1868A	Neonatal Single Patient NIBP Cuff, Size 2	55
M1870A	Neonatal Single Patient NIBP Cuff, Size 3	55

Philips Reusable SpO₂ Sensors

M1193A	Reusable Neonatal Wrap Sensor	61
M1193T	Reusable Neonatal Wrap Sensor	61
M1195A	Reusable Infant Finger Glove Sensor	62

SpO₂ Adapter Cables

989803148221	Masimo® MP10 Adapter Cable (LNCS sensors)	63
M1020-61102	Masimo LNOP Adapter Cable (12-pin)	63
M1943NL	OxiMax® SpO ₂ Adapter Cable (8-pin)	64

Product removed**Multi Patient NIBP Cuff Kits – Single Hose**

864305	Neo IEC Reusable Starter Kit – New NIBP	19
864306	Neo AAMI Reusable Starter Kit – New NIBP	19
864307	Neo IEC Disposable Starter Kit – New NIBP	19
864308	Neo AAMI Disposable Started Kit – New NIBP	19

Dear Customer:

If you're looking for high-quality medical supplies to optimize the performance of your medical devices, you'll find them right here. Whether your focus is on quality patient care, reliability, or cost containment, you can turn to a single source – Philips Medical Supplies. Explore our complete family of supplies, from electrodes to SpO₂ sensors, for every area of care.

With Philips, you can buy your medical equipment and supplies from a single source, simply and easily. Use this guide to find the Philips supplies you need – or go to our website to purchase supplies: <http://www.philips.com/healthcarestore>.

Using this guide

Latex Cautions

Philips Medical Supplies currently distributes over 600 products, of which less than 9% contain natural rubber latex or dry natural rubber that could contact humans. Because these substances may cause allergic reactions, products in this guide that contain them are identified with the following symbol:

LATEX

Caution: Products that are identified with this symbol contain dry natural rubber or natural rubber latex that can come into contact with humans and may cause allergic reactions.

To the best of our knowledge, products not identified with this symbol are not known to contain natural rubber latex or dry natural rubber. The latex symbol indicates **only** which **medical supplies** contain natural rubber latex or dry natural rubber that can come into contact with humans and may cause allergic reactions. If you want to know whether a Philips **medical device** contains natural rubber latex or dry natural rubber, refer to the product labeling.

The FDA requires that products be labeled with human contacting latex content information when appropriate. Always refer to product labeling for accurate information on human contacting dry natural rubber or natural rubber latex.

CE Marking

Products labeled with this symbol conform to the European Union (EU) directives as mandated by EU member state legislation and are available for sale in Europe to countries that accept CE marking with restrictions as may be noted. Medical devices labeled with this marking means the

product is in conformity with the European Medical Directive (MDD): 93/42/EEC Council Directive of 14 June 1993 on the approximation of the laws of the Member States concerning medical devices.

In order for a product to be able to be labeled with the CE marking, product quality and manufacturing quality systems must be documented to comply with certain standards. Compliance with the European directives and regulations ensures a standard level of performance for device function, safety, and health protection for users, patients, and others.

The CE symbol provided with some products in the Philips Medical Supplies Reference Guide and online store is intended only as an indication of availability in EU countries and is not intended as a legal product marking.

Countries that accept CE marked products include Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Russia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom. A CE mark does not necessarily indicate that a product is marketed or distributed in all of the countries listed above.

Non-medical products that comply with other harmonized directives may also bear the CE mark.

Note: Some products identified in this catalogue may not be licensed in accordance with Canadian law.

Ordering Information

To order the medical supplies in this guide or to obtain more information or additional literature, contact your Philips Medical Supplies representative near you or visit our website:

<http://www.philips.com/healthcarestore>

Online Ordering

Online ordering is a simple process:

1. Log in using a secure password.
2. Enter the product number and quantity in the shopping cart.
3. Enter your PO or credit card number at check-out.
4. Your order is immediately transferred into our fulfillment system, where you can track it anytime.

If you are new to the web, we can guide you through the process on the phone in minutes. Many new users tell us that the ordering process is so intuitive that they do not need help, even the first time.

You can search for supplies in the online catalog three easy ways: by application category, Philips instrument, or keyword.

If you currently transmit via MMS/EDI or auto fax, the interactive system provides real-time information and services not available through one-way transmissions.

Online ordering is available to U.S. customers at the following url:

<http://www.philips.com/newhealthcarestore>

OB/Neonatal Monitoring Supplies

Fetal Monitoring Supplies

Fetal Belts	9
Fetal Direct ECG Monitoring	12
Fetal Paper	14
Intrauterine Pressure (IUP) Monitoring	15
Miscellaneous Supplies	16

Maternal Monitoring Supplies

Non-Invasive Blood Pressure	19
Pulse Oximetry	26

Neonatal Monitoring Supplies

Capnography	29
Spirometry	34
Transcutaneous Gas Monitoring	35
ECG Cables and Leads	37
ECG Electrodes	42
Invasive Blood Pressure	45
Non-Invasive Blood Pressure	47
Pulse Oximetry	55
Temperature	61



Fetal Monitoring Supplies

Fetal Belts

Belt Guides and Buttons

Product #	Description	Packaging	Use With:
989803143401	Belt Clip for TOCO+ Plastic Belt Clip for TOCO+ Transducer	6 clips	Toco+ Transducer M2735A; Toco Transducer M2734A, M2734B; US Transducer M2736A
15155-43101	Abdominal Belt Guide	1 guide	15245A US Transducer; 15248A Toco Transducer
M1569A	Fetal Belt Buttons	10 buttons	M2734B Toco Transducer, M2734A, M2735A, M2736A

Disposable Belts

Product #	Description	Packaging	Use With:
14328A CE	Disposable Abdominal Belt: white foam with velcro, 1.4m long, 50mm wide	50 belts	15245A US Transducer; 15248A Toco Transducer
M2208A CE	Disposable Abdominal Belt: yellow, with buttons, 1.3m long, 60mm wide	100 rolls	M1356A, US Transducer; M1355A, M2734A, M2734B, M2735A, M2736A Toco/Toco+ Transducer; M1356A, M2736A US Transducer; M1358A MECG/US Transducer; M1364-60001 ECG Module; M1365A FSpO ₂ /ECG Module; M2725A Wireless Toco Transducer; M2726A Wireless US Transducer; M2727A Wireless ECG Transducer; M2738A ECG/IUP Module, M2720A



989803143401



14328A



M2208A

Reusable Belts

Product #	Description	Packaging	Use With:
1500-0628 LATEX	Reusable Leg Belt: brown, washable, bulk roll, contains Latex, 15m long, 32mm wide	1 bulk roll	M1364-6001 ECG Module, M1365A FSpO ₂ /ECG Module, M1357A DECG Module, M1359A MECG Module, 15240A DECG Transducer, M2738A ECG/IUP Module
1500-0642 LATEX	Reusable Abdominal Belt: brown, washable, pre-cut, contains Latex, 1.3m long, 60mm wide	5 rolls	M1356A, M2736A US Transducer; M2734A, M2734B, M2735A Toco/Toco+ Transducer; M1358A MECG/US Transducer; M1364-60001 ECG Module; M1365A FSpO ₂ /ECG Module; M2738A ECG/IUP Module; M2725A Wireless Toco Transducer; M2726A Wireless US Transducer; M2727A Wireless ECG Transducer, M2720A
1500-0643 LATEX	Reusable Abdominal Belt: brown, washable, bulk roll, contains Latex, 15m long; 60mm wide	1 bulk roll	M1356A, M2736A US Transducer; M1355A, M2734A, M2734B, M2735A Toco/Toco+ Transducer; M1358A MECG/US Transducer; M1364-60001 ECG Module; M1365A FSpO ₂ /ECG Module; M2738A ECG/IUP Module; M2725A Wireless Toco Transducer; M2726A Wireless US Transducer; M2727A Wireless ECG Transducer, M2720A
M1562A CE LATEX	Reusable Abdominal Belt: brown, washable, pre-cut, buttons attached, contains Latex, 1.3m long, 50mm wide	5 rolls	15245A, M1356A, M2726A, M2736A US Transducer; 15248A, M1355A, M2734A, M2734B, M2735A Toco/Toco+ Transducer; M1357A DECG Transducer; M1358A, M1359A MECG/US Transducer; M1364-60001 ECG Module; M1365A FSpO ₂ /ECG Module; M2738A ECG/IUP Module; M2725A Wireless Toco Transducer; M2726A Wireless US Transducer; M2727A Wireless ECG Transducer, M2720A



1500-0642



M1562A

Reusable Belts (continued)

Product #	Description	Packaging	Use With:
M1562B CE	Reusable Abdominal Belt: grey, washable, water resistant, pre-cut, 1.3m long, 50mm wide	5 rolls	15245A, M1356A, M2726A, M2736A US Transducer; 15248A, M1355A, M2734A, M2734B, M2735A Toco/Toco+ Transducer; M1357A DECG Transducer; M1358A, M1359A MECG/US Transducer; M1364-60001 ECG Module; M1365A FSpO ₂ /ECG Module; M2738A ECG/IUP Module; M2725A Wireless Toco Transducer; M2726A Wireless US Transducer; M2727A Wireless ECG Transducer, M2720A
M4601A CE	Reusable Leg Belt: grey, washable, water resistant, bulk roll, 15m long, 32mm wide	1 bulk roll	M1364-6001 ECG Module, M1365A FSpO ₂ /ECG Module, M1357A DECG Module, M1359A MECG Module, 15240A DECG Transducer, M2738A ECG/IUP Module
M4602A CE	Reusable Abdominal Belt: grey, washable, water resistant, pre-cut, 1.3m long, 60mm wide	5 rolls	15245A, M1356A, M2736A US Transducer; 15248A, M1455A, M2734A, M2734B, M2735A Toco/Toco+ Transducer; M1358A MECG/US Transducer; M1364-60001 ECG Module; M1365A FSpO ₂ /ECG Module; M2738A ECG/IUP Module; M2725A Wireless Toco Transducer; M2726A Wireless US Transducer; M2727A Wireless ECG Transducer, M2720A
M4603A CE	Reusable Abdominal Belt: grey, washable, water resistant, bulk roll, 15m long, 60mm wide	1 bulk roll	



M1562B



M4601A



M4603A

Fetal Direct ECG Monitoring

Philips DECG Legplate Adapter Cables

Product #	Description	Packaging	Use With:
989803137651 CE	DECG Reusable Legplate Adapter Cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	Fetal Monitors: M2703A, M2705A, M1350A, M1350B, M1350C, M1353A, Transducers/ Modules: M1357A, M1364-60001, M2738A, M1365A, M2727A, M2735A, M2720A
989803137791 CE	DECG Reusable Legplate Adapter Cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	GE Corometrics Monitors: 115, 116, 118, 119, 173, 174; GE Series Monitors :120, 126, 128, 129, 250 using a GE Corometrics Y-adapter cable 1442AAO or socket adapter 1442BAO.
989803137801 CE	DECG Reusable Legplate Adapter Cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	GE Corometrics Monitors: 120, 126, 128, 129
989803137811 CE	DECG Reusable Legplate Adapter Cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	Spacelabs IM77 (other models of Spacelabs may require an additional adapter cable.)
989803139771 CE	Attachment Electrode, minimum shelf life 6 months	25 electrodes	Fetal Monitors: M2705A, M2703A, M1350A, M1350B, M1350C, M1353A, Transducers/ Modules: M1357A, M1364-60001, M2738A, M1365A, M2727A, M2735A, M2720A



989803137651



989803137791



989803137801




989803137811



989803139771

Philips Fetal Scalp Electrodes

Product #	Description	Packaging	Use With:
989803137631 CE	Fetal Spiral Electrode, single spiral	25 electrodes	Use with DECG Legplate Adapter Cables: 989803137651, 989803137791, 989803137801, 989803137811, M2735A, M2738A, M2720A, M1353A, M1350B, M2703A, M2705A
989803137641 CE	Fetal Spiral Electrode, double spiral (available in Europe only)	25 electrodes	Use with DECG Legplate Adapter Cables: 989803137651, 989803137791, 989803137801, 989803137811, M2735A, M2738A, M2720A, M1353A, M1350B, M2703A, M2705A



989803137631




989803137641

Fetal Paper

Chemical/Thermal Paper

Product #	Description	Packaging	Use With:
9270-0485 CE	Fetal Recording Paper: orange grid, 1-min. (3cm) highlights, 1-100 mmHg IUP scale, 30-240 FHR scale.	40 packs (150 pages per pack)	Fetal Monitors: 8030A, 8031A, 8032A, 8032B, 8040A, 8041A
9270-0630-050 CE	Fetal Recording Paper: green grid, 1-100 mmHg IUP scale, 0-13 kPa scale, 50-210 FHR scale	50 packs (150 pages per pack)	Fetal Monitors: 8030A, 8031A, 8040A, 8041A
M1910A CE	Fetal Recording Paper: orange grid, 1-min. (3cm) highlights, 0-100 mmHg IUP scale, 30-240 FHR scale.	40 packs (150 pages per pack)	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350A, M1350B, M1350C, M1351A, M1353A, Recorder M2925A
M1911A CE	Fetal Recording Paper: green grid, 0-100 mmHg IUP scale, 0-13 kPa scale, 50-210 FHR scale	40 packs (150 pages per pack)	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350A, M1350B, M1350C, M1351A, M1353A, Recorder M2925A
M1913A CE	Fetal Recording Paper: green grid, 1-min. (3cm) highlights, 0-100 mmHg IUP scale, 50-210 FHR scale.	40 packs (150 pages per pack)	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350A, M1350B, M1350C, M1351A, M1353A, Recorder M2925A



M1910A



M1911A

Intrauterine Pressure (IUP) Monitoring

Disposable Koala™ IUP Solution

Product #	Description	Packaging	Use With:
989803143931 CE	Reusable IUP Connector Cable (Avalon) for Koala™ disposable IUP Catheter	1 cable	M2738A Patient Module, M2735A TOCO+ Transducer, Fetal Monitors M2703A, M2705A, M1350A/B/C
M1333A CE	Koala™ disposable IUP Catheters	10 catheters	Fetal Monitors: M2703A, M2705A, 8040A, M1350A, M1350B, M1350C, M1351A, M1353A, M2735A, M2738A
M1334A CE	Reusable IUP Connector Cable for Koala™ disposable IUP Catheter	1 cable	Fetal Monitors: 8040A, M1350A, M1350B, M1350C, M1351A, M1353A

Reusable Memscap™ IUP Solution

Product #	Description	Packaging	Use With:
CPJ84022 CE	Single Use Sterile Domes	50 domes	CPJ840J5, CPJ840J6, M1310A, M1350B, M1350C, M1353A, M3015A
CPJ84026 CE	Monitoring Line Set	10 sets per bag	863065, 863066, 863068, 863087, 863088, M3536M, M3002A, M8102A, M3535A, M8105AT, M8105A, M3536A, M3015A, M3015B
CPJ84046 CE	Transducer holder	4 holders	CPJ840J5, M3015A
CPJ840J5 CE	Memscap Reusable Pressure Transducer for IUP monitoring, Cable length: 3 meters, includes one transducer holder.	1 transducer and holder	M1310A, M1350B, M1350C, M1353A
CPJ84447 CE	Pressure Transducer Pole Mount		CPJ840J5 transducer with instruments: M1350B, M1350C, M1353A, M3015A



989803143931



M1333A



M1334A

Miscellaneous Supplies

Miscellaneous Fetal Supplies

Product #	Description	Packaging	Use With:
989803143411 CE	Avalon FM30/50-Remote Event Marker	1 event marker	M2702A, M2703A, M2704A, M2705A Avalon fetal monitors
M4503A	Fetal Record File Kit: 5 sets of tabs, 60 record wrappers, 120 patient labels, 40 antepartum labels, 60 delivery labels, 1 file box label	1 kit	

Modules, Transducers, Adapters

Product #	Description	Packaging	Use With:
M2738A CE	Avalon Patient Module: for IUP, fetal/maternal ECG	1 module	Avalon Fetal Monitors: M2702A, M2703A, M2704A, M2705A



989803143411



M2738A

Fetal Supplies Kits

Product #	Description	Packaging	Use With:
862214	Standard Antepartum Kit (available in N. America only)	(10) M1569A (40) M1910A (12) 40483A gel (100) M2208A belts (6) Toco+ clips Avalon Fetal	Avalon Fetal Monitors, FM20, FM30, FM40, FM50, M2702A, M2703A, M2704A, M2705A
862215 CE	Standard Antepartum Kit (available in CN, EU, LA, Asia (except Japan) only)	(10) M1569A (40) M1911A (12) 40483A gel (100) M4602A belts (6) Toco+ clips	Avalon Fetal Monitors, FM20, FM30, FM40, FM50, M2702A, M2703A, M2704A, M2705A
862217	Intrapartum DECG Kit (Single Spiral) (not available in CN, MX, RU)	(25) 989803137631 FSE (25) 989803139771 electrodes	Avalon Fetal Monitors, FM30, FM50, M2703A, M2705A, M1353A, M1350B
862218 CE	Intrapartum DECG Kit (Double Spiral) (available in EU only)	(25) 989803137641 FSE (25) 989803139771 electrodes	Avalon Fetal Monitors, FM30, FM50, M2703A, M2705A, M1353A, M1350B

Disposable ECG Electrodes

Product #	Description	Packaging	Use With:
40493D CE	Adult foam ECG electrode	5 per pack 300 per case	Fetal/Maternal Monitors: 8040A, M1350B, M1350C, M2702A, M2703A, M2705A, M2704A, M2720A
40493E CE	Adult foam ECG electrode	30 per pack (300 per case)	Fetal/Maternal Monitors: 8040A, M1350B, M1350C, M2702A, M2703A, M2705A, M2704A, M2720A

ECG Lead Set

Product #	Description	Packaging	Use With:
M1363A CE	MECG cable with reusable leadset for maternal ECG	1 cable	Fetal Monitors: Patient Module: M1364A, M1365A, M2727A, M2738A, M2735A, M2702A, M2703A, M2704A, M2705A, M2720A

Ultrasound Transmission Gel

Product #	Description	Packaging	Use With:
40483A CE	Ultrasound Transmission Gel, 8 1/2oz refillable squeeze bottle Shelf life: at least 6 months	12 bottles	Ultrasound transducers, M2720A Avalon Wireless Transducer System
40483B CE	Ultrasound Transmission Gel, 5 liter bottle use to refill 40483A. Shelf life: at least 6 months	(1) 5 liter bottle	Ultrasound transducers, M2720A Avalon Wireless Transducer System



40493D-E



M1363A



40483A

Maternal Monitoring Supplies

Non-Invasive Blood Pressure

Reusable NIBP Cuffs

Multi-Patient Comfort Cuffs – Single Hose

Product #	Description	Packaging	Use With:
REVISED M1573A CE	Small Adult, light blue tubing Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference = 20.5cm to 28cm	1 cuff	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M1574A CE	Adult, navy blue tubing Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27cm to 35cm	1 cuff	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M1575A CE	Large Adult, burgundy tubing Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference = 34cm to 43cm	1 cuff	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M1576A CE	Thigh, grey tubing Bladder width: 7.9" (20cm) Bladder length 16.9" (43cm) Limb circumference = 42cm to 54cm	1 cuff	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



M1573A – M1576A

Multi-Patient Easy Care Cuffs – Single Hose

Product #	Description	Packaging	Use With:
REVISÉ M4554B/ M4554B5 CE	Small Adult, light blue. Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5cm to 28.5cm	1 cuff/5 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ M4555B/ M4555B5 CE	Adult, navy blue. Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5cm to 36cm	1 cuff/5 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ M4556B/ M4556B5 CE	Adult Extra-Long, navy blue. Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5cm to 36cm	1 cuff/5 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ M4557B/ M4557B5 CE	Large Adult, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35cm to 45cm	1 cuff/5 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ M4558B/ M4558B5 CE	Large Adult Extra-Long, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35cm to 45cm	1 cuff/5 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ M4559B/ M4559B5 CE	Thigh, grey. Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference = 44cm to 56cm	1 cuff/5 cuffs	Fetal Monitors: M1350B, M1350C, M2702A, M2703A, M2704A, M2705A, M8105AS, 863275, 863276, 863278



M4552B – M4559B

Multi-Patient NIBP Cuff Kits – Single Hose

Product #	Description	Packaging	Use With:
REVISÉ 864288 CE	Easy Care Adult Kit – 4 sizes	consists of 4 cuffs: M4554B, M4555B, M4557B, M4559B	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ 864290 CE	Easy Care Assortment Kit – 6 sizes	consists of 6 cuffs: M4552B, M4553B, M4554B, M4555B, M4557B, M4559B	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ 864291 CE	Easy Care Adult Long Kit – 6 sizes	consists of 6 cuffs: M4554B, M4555B, M4556B, M4557B, M4558B, M4559B	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISÉ M1578A CE	Multi-Patient Comfort Cuff Assortment Kit	consists of 4 cuffs: M1573A, M1574A, M1575A, M1576A	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A

Disposable NIBP Cuffs

Multi Patient Disposable Cuffs – Multi Care

Product #	Description	Packaging	Use With:
REVISED 989803183331 CE	Small Adult, light blue Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference = 20.5cm to 28cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED 989803183341 CE	Adult, navy blue Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27cm to 35cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
NEW 989803183351 CE	Adult XL (navy blue with edge striping) Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27cm to 35cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED 989803183361 CE	Large Adult, burgundy Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference = 34cm to 43cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED 989803183371 CE	Thigh, grey Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm) Limb circumference: 42cm to 54cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



Gentle Care Disposable Cuffs

Product #	Description	Packaging	Use With:
REVISED M4574B CE	Small Adult, light blue. Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5cm to 28.5cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M4575B CE	Adult, navy blue. Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5cm to 36cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M4576B CE	Adult Extra-Long, navy blue. Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5cm to 36cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M4577B CE	Large Adult, burgundy. Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35cm to 45cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M4578B CE	Large Adult Extra-Long, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35cm to 45cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M4579B CE	Thigh, grey. Bladder width: 8.3" (21cm) Bladder length 18.9" (48cm) Limb circumference = 44cm to 56cm	10 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



M4574B – M4579B

Single Care Disposable Cuffs

Product #	Description	Packaging	Use With:
NEW 989803182291 CE	Small Adult (light blue) Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5cm to 28.5cm	20 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
NEW 989803182301 CE	Adult (navy blue) Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5cm to 36cm	20 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
NEW 989803182311 CE	Adult XL (navy blue with edge striping) Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5cm to 36cm	20 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
NEW 989803182321 CE	Large Adult (burgundy) Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35cm to 45cm	20 cuffs	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



NIBP Air Hoses and Caps

NIBP Cuff Caps

Product #	Description	Packaging	Use With:
989803151601 CE	Cuff Caps for Single Tube Cuff (10) NIBP cuff caps, single hose	10	n/a

Adult NIBP Pressure Air Hose

Product #	Description	Packaging	Use With:
REVISED M1598B CE	Adult pressure air hose 1.5m length: 4.92' (1.5m)	1 air hose	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
REVISED M1599B CE	Adult pressure air hose 3.0m length: 9.84' (3.0m)	1 air hose	Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



989803151601



M1598B



M1599B

Pulse Oximetry

Reusable SpO₂ Sensors

Philips Reusable SpO₂ Sensors

Product #	Description	Packaging	Use With:
REVISED M1191B CE	Adult Finger Glove Sensor Length = 2m Patient size: >110lb (>50kg) 8 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. To add length, use optional SpO ₂ Extension Cable: M1941A Fetal Monitors: M1350B, M1350C Requires the SpO ₂ Adapter Cable: M1940A
REVISED M1191BL CE	Adult Finger Glove Sensor Length = 3m Patient size: >110lb (>50kg) 8 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A No cable required. M8105AS, 863276, 863278 Extension cable not allowed
REVISED M1191T CE	Adult Finger Glove Sensor Length = 45cm Patient size: >110lb (>50kg) 9 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. Requires SpO ₂ Adapter Cable: M1943A/AL Fetal Monitors: M1350B, M1350C Requires SpO ₂ Adapter Cable: M1900B
REVISED M1192A CE	Pediatric/Small-Adult Glove Sensor Length = 1.5m Patient size: 33 – 110lb (15 – 50kg) 8 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. To add length, use optional SpO ₂ Extension Cable: M1941A Fetal Monitors: M1350B, M1350C Requires the SpO ₂ Adapter Cable: M1940A
REVISED M1192T CE	Pediatric/Small-Adult Glove Sensor Length = 45cm Patient size: 33 – 110lb (15 – 50kg). 9 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. Requires SpO ₂ Adapter Cable: M1943A/AL



M1191B/BL



M1191T



M1192A/M1192T

Philips Reusable SpO₂ Sensors (continued)

Product #	Description	Packaging	Use With:
REVISED M1194A CE	Pediatric/Adult Ear Clip Sensor Length = 1.5m Patient size: >88lb (>40kg) 8 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. To add length, use optional SpO ₂ Extension Cable: M1941A Fetal Monitors: M1350B, M1350C Requires the SpO ₂ Adapter Cable: M1940A
REVISED M1196A CE	Pediatric/Adult Finger Clip Sensor Length 3m Patient size: >88lb (>40kg) 8 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 860335, 860338, 863276, 863278. Extension cable not allowed
REVISED M1196S CE	Reusable Clip Adult SpO ₂ Sensor (2m) 2m (79") cable, Philips connector, finger Patient size: >40kg (>88lb)	1 sensor	Fetal Monitors: M2703A, M2704A, M2705A, M8102A, M8105A, M8105AT, M8105AS, 863276, 863278
REVISED M1196T CE	Pediatric/Adult Finger Clip Sensor Length 90cm Patient size: >88lb (>40kg) 9 pin connector	1 sensor	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. Requires SpO ₂ Adapter Cable: M1943A/AL Fetal Monitors: M1350B, M1350C Requires SpO ₂ Adapter Cable: M1900B
REVISED M1627A CE	Wristband Velcro straps for SpO ₂ glove sensors	10 straps	



M1194A






M1196A



M1196T

Disposable SpO₂ Sensors

Philips Disposable SpO₂ Sensors

Product #	Description	Packaging	Use With:
REVISED M1131A 	Adult/Pediatric Finger Clip Disposable Sensor (single patient use) Length = 45cm. Patient size: >44lb (>20kg) 9-pin (D-sub) connector	20 sensors	Fetal Monitors: M2703A, M2705A Requires SpO ₂ Adapter Cable: M1943A/AL, M8105AS, 860335, 860338, M1350B, M1350C, 863276, 863278
REVISED M1133A 	Neonatal/Infant/Adult Disposable Wrap Sensor (single patient use) Length = 90cm Patient sizes: adult >88lb (>40kg) 9-pin (D-sub) connector	20 sensors	Fetal Monitors: M2703A, M2705A Requires SpO ₂ Adapter Cable: M1943A/AL, M8105AS, 860335, 860338, M1350B, M1350C, 863276, 863278
REVISED M1134A 	Adhesive-Free Neonatal/Infant/Adult Disposable Wrap Sensor (single patient use) Length = 90cm Patient size: adult >88lb (>40kg) 9-pin (D-sub) connector	20 sensors	Fetal Monitors: M2703A, M2705A Requires SpO ₂ Adapter Cable: M1943A/AL, M8105AS, 860335, 860338, M1350B, M1350C, 863276, 863278



M1131A




M1133A



M1134A

Adhesive SpO₂ Sensors

Product #	Description	Packaging	Use With:
REVISED M1904B 	Adult Disposable Sensor Patient size: greater than 66lb (>30kg) 9-pin (D-sub) connector (not available in N. America and Japan)	24 sensors and tapes	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278. Requires SpO ₂ Adapter Cable: M1943A/AL Fetal Monitors: M1350B, M1350C Requires SpO ₂ Adapter Cable: M1900B.



M1904B

SpO₂ Sensor Adapter Cables

Philips SpO₂ Adapter Cables

Product #	Description	Packaging	Use With:
M1900B CE	SpO ₂ Adapter Cable (12-pin) Length = 9.8' (3m) Adapts 9-pin sensors to 12-pin sockets	1 cable	Fetal Monitors: M1350B, M1350C
M1940A CE	SpO ₂ Adapter Cable (12-pin) Length = 6.6' (2m) Adapts 8-pin sensors to 12-pin sockets	1 cable	Fetal Monitors: M1350B, M1350C
REVISED M1941A CE	SpO ₂ Extension Cable, 6.6' (2m) For 8-pin sensors to 8-pin sockets Adds additional length to sensor cable	1 cable	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278
M1943A CE	SpO ₂ Adapter Cable (8-pin) Length = 3.3' (1.1m) Adapts 9-pin sensors to 8-pin sockets	1 cable	Fetal Monitors: M2703A, M2705A, M8105AS, 860338, 860335, 863276, 863278
M1943AL CE	SpO ₂ Adapter Cable (8-pin) Length = 9.8' (3m) Adapts 9-pin sensors to 8-pin sockets	1 cable	Fetal Monitors: M2703A, M2705A, M8105AS, 863276, 863278



M1900B



M1940A



M1941A



M1943A/M1943AL

Neonatal Monitoring Supplies

Capnography

Sidestream CO₂

Microstream® FilterLine®

Product #	Description	Packaging	Use With:
989803159581 CE	VitaLine H Set Infant/Neonatal, Intubated, long term. High ambient humidity environments. Patient size: (<2kg)	25 sample sets	Microstream Ext. Server: M3015A Monitors: VM1 863266, VM8 863066, MP5 M8105A, M8105AS, M3015B
989803160261 CE	FilterLine H Set (long) Infant/Neonatal, Intubated, short term. Length = 4m, Patient size: (<2kg)	25 sample sets	Microstream Ext. Server: M3015A Monitors: VM1 863266, VM8 863066, MP5 M8105A, M8105AS, M3015B
989803159571 CE	VitaLine™ H Set, Adult/Pediatric	25 sets	Microstream Ext: M3015A, VM1: 863266; MP5: M8105A, M8105AS, M3015B
989803159581 CE	VitaLine™ H Set, Infant/Neonatal	25 sets	Microstream Ext: M3015A, VM1: 863266; MP5: M8105A, M8105AS, M3015B
989803160241 CE	FilterLine® Set Long, Adult/Pediatric	25 sets	Microstream Ext: M3015A, VM1: 863266; MP5: M8105A, M8105AS, M3015B
989803160251 CE	FilterLine® H Set Long, Adult/Pediatric	25 sets	Microstream Ext: M3015A, VM1: 863266; MP5: M8105A, M8105AS, M3015B
989803160261 CE	FilterLine® H Set Long, Infant/Neonatal	25 sets	Microstream Ext: M3015A, VM1: 863266; MP5: M8105A, M8105AS, M3015B
989803160271 CE	Smart CapnoLine® O ₂ Long, Pediatric	25 cannula	
989803177971 CE	Smart CapnoLine® H O ₂ Pediatric	25 Filterlines	
989803177981 CE	Smart CapnoLine® H O ₂ Pediatric Long	25 Filterlines	
989803178001 CE	CapnoLine® H O ₂ Infant/Neonatal	25 Filterlines	
989803178011 CE	CapnoLine® H Infant/Neonatal Long	25 Filterlines	



989803159581



989803160261

Microstream® FilterLine® (continued)

Product #	Description	Packaging	Use With:
989803178021 CE	Nasal FilterLine® Infant/Neonatal	25 Filterlines	
989803178031 CE	Smart CapnoLine® Guard	25 Filterlines	
989803178041 CE	Smart CapnoLine® Guard O ₂	25 Filterlines	
989803178051 CE	Smart CapnoLine® Guard O ₂ Long	25 Filterlines	M3015A, M3015B, M8105A, M8105AS
989803178071 CE	Hook and Loop Strap	1 box = 25	
989803179121 CE	Nasal Filterline® O ₂ Pediatric	25 Filterlines	M3015A, M3015B, M8105A, M8105AS
M1920A CE	FilterLine® Set, Adult/Pedi, Intubated	25 sets	Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A, M8105AS, M3015B
M1921A CE	FilterLine® H Set, Adult/Pedi, Intubated	25 sets	Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A, M8105AS, M3015B
M2520A CE	Smart CapnoLine® O ₂ , Pediatric, Non-intubated Dual Purpose	25 cannula	Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A, M8105AS, M3015B
M2524A CE	Smart CapnoLine®, Pediatric, Non-intubated Single Purpose	25 cannula	Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A, M8105AS, M3015B
M4681A CE	CapnoLine® H O ₂ , Pediatric, nasal, (up to 24 hours) Non-intubated Single Purpose	25 cannula	Microstream Ext: M3015A, VM1: 863266, MP5: M8105A, M8105AS, M3015B
M4687A CE	NIV Line, Pediatric, nasal (up to 8 hours), mask, single purpose	25 lines	Microstream Ext: M3015A, VM1: 863266, MP5: M8105A, M8105AS, M3015B
M1923A CE	FilterLine H Set Infant/Neonatal, Intubated, long term. Patient size: (<2kg)	25 sample sets	Microstream Ext. Server: M3015A Monitors: VM1 863266, VM8 863066, MP5 M8105A, M8105AS, M3015B
M4691A CE	Capnoline H, Infant/Neonatal, Non-intubated, single purpose, long term (nasal). Patient size: (<10kg)	25 cannula	Microstream Ext. Server: M3015A Monitors: VM1 863266, VM8 863066, MP5 M8105A, M8105AS, M3015B



M1923A



M4691A

Philips IntelliVue Sidestream CO₂

Product #	Description	Packaging	Use With:
M2741A CE	Sidestream CO ₂ Sensor	1 sensor	M3014A Capnography Extension, Monitors: M3002A X2, M8102A MP2, M8105A MP5, 860335, 860338
M2745A CE	CO ₂ Nasal Cannula - Pediatric, disposable, non-intubated	10 cannula	Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A, 860335, 860338
M2746A CE	Single Patient Disposable CO ₂ Nasal Cannula, non-intubated. Patient size: Infant	10 cannula	M2741A Sidestream Sensor with Sidestream Ext. Server: M3014A. Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A, 860335, 860338
M2751A CE	CO ₂ / O ₂ Nasal Cannula - Pediatric, disposable, non-intubated	10 cannula	Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A, 860335, 860338
989803144471 CE	Single Patient Disposable CO ₂ / O ₂ Nasal Cannula, non-intubated. Patient size: Infant	10 cannula	M2741A Sidestream Sensor with Sidestream Ext. Server: M3014A. Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A
M2757A CE	CO ₂ Oral-Nasal Cannula - Pediatric, disposable, non-intubated	10 cannula	Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A, 860335, 860338



M2741A



M2746A



989803144471

Philips IntelliVue Sidestream CO₂ (continued)

Product #	Description	Packaging	Use With:
M2761A CE	CO ₂ / O ₂ Oral-Nasal Cannula - Pediatric, disposable, Non-intubated	10 cannula	Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A, 860335, 860338
M2776A CE	Straight Sample Line, Intubated, non-humidified. Recommended for HME filter.	10 lines	M2741A Sidestream Sensor with Sidestream Ext. Server: M3014A. Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A
M2777A CE	Straight Sample Line H, Intubated, humidified. Recommended for HME filter.	10 lines	M2741A Sidestream Sensor with Sidestream Ext. Server: M3014A. Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A
M2773A CE	Airway Adapter Set H, Intubated (ET equal to or less than 4.0mm)	10 airway adapters	M2741A Sidestream Sensor with Sidestream Ext. Server: M3014A. Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A, 860335, 860338
989803144531 CE	Airway Adapter Set - ET <=4.0mm	10 airway adapters	M2741A Sidestream Sensor with Sidestream Ext. Server: M3014A. Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A, 860335, 860338



M2776A



M2777A



M2773A



989803144531

Mainstream CO₂

Philips IntelliVue Mainstream CO₂ Monitoring

Product #	Description	Packaging	Use With:
M2501A CE	Mainstream CO ₂ Sensor	1 sensor	M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A, 860335, 860338
M2513A CE	Reusable airway adapter, Adult/Pediatric (color: black) ET tube size: greater than 4mm, deadspace less than 5cc	1 adapter	M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A, 860335, 860338
M2516A CE	Reusable CO ₂ airway adapter (color: red) ET tube size: less than 4mm, deadspace less than 1cc	1 airway adapter	M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A
M2533A CE	Single patient disposable CO ₂ airway adapter, Adult/Pediatric (color: clear) ET tube size: greater than 4mm, deadspace less than 5cc	10 adapters	M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A, 860335, 860338
M2536A CE	Single Patient Disposable CO ₂ airway adapter (color: violet) ET tube size: less than 4mm, deadspace less than 1cc	10 airway adapters	M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, MMS X2: M3002A, 860335, 860338



M2501A



M2516A



M2536A

Spirometry

Flow Sensors and Combined CO₂/Flow

Product #	Description	Packaging	Use With:
M2781A CE	Combined CO ₂ / Flow Sensor. Adult / Pediatric, disposable ET: greater than 5.5mm (color: clear)	10 sensors	Spirometry Module: M1014A with M2501A Mainstream CO ₂ Sensor and M3014A Capnography Ext.
M2782A CE	Combined CO ₂ / Flow Sensor Neonatal. Infant, disposable ET: 2.5mm to 4mm (color: violet)	10 sensors	Spirometry Module: M1014A with M2501A Mainstream CO ₂ Sensor and M3014A Capnography Ext.
M2783A CE	CO ₂ / Flow Sensor Pediatric, disposable ET: 3.5mm to 6mm (color: green)	10 sensors	Spirometry Module: M1014A with M2501A Mainstream CO ₂ Sensor and M3014A Capnography Ext.
M2785A CE	Flow Sensor Adult/Pediatric, disposable ET: greater than 4mm (color: clear)	10 sensors	M1014A Spirometry Module
M2786A CE	Flow Sensor, infant/neonatal ET: less than 4mm. (color: violet)	10 flow sensors	M1014A Spirometry Module



M2782A



M2786A

Transcutaneous Gas Monitoring

tcpO₂/tcpCO₂

Product #	Description	Packaging	Use With:
M1918A CE	Solid state tcpO ₂ /tcpCO ₂ transducer	Includes: 1 transducer, 1 cap	Module: M1018A
M1918B CE	Solid state tcpO ₂ /tcpCO ₂ transducer	Includes: 1 transducer, 1 Cap, 1 Calibration Ring	Module: M1018A
989803151951	tcpO ₂ /tcpCO ₂ Membrane Kit	10 tools containing membranes, 10ml electrolyte, 2 sensor caps, 5 cleaning papers	M1918B transducer
989803151961	tcpO ₂ /tcpCO ₂ Fixation Kit	100 fixation rings, 4 bottles contact fluid @ 20ml	M1918B transducer
989803160371	Calibration Ring	1 ring	M1018A
15209-60010 CE	Accessory Kit for tcpO ₂ /tcpCO ₂ Transducer	10ml electrolyte, 12 membranes, 2 transducer caps, 25 cleaning papers, Instructions, O-ring remover	M1918A transducer
15209-60020 CE	Fixation kit for tcpO ₂ /tcpCO ₂ transducers	100 fixation rings, 4 (20ml) bottles contact liquid, instructions	M1918A transducer



M1918B



989803151951



989803151961



989803160371

tcpO₂/tcpCO₂ (continued)

Product #	Description	Packaging	Use With:
15210-60010	Gas Cal 1 Gas Cylinders for tcpO ₂ /tcpCO ₂ calibration. (not available in Germany)	6 cylinders	15205A, 15206A, 15209A, 15209AH, 15210A, 15210B, 15210AH, 15210B, M1018A, M1918A, M1918B
15210-60020	Gas CAL 2 gas cylinders for tcpO ₂ /tcpCO ₂ calibration. (not available in Germany)	6 cylinders	15205A, 15206A, 15209A, 15209AH, 15210A, 15210B, 15210AH, M1018A
15210-64010	Gas Cal 1 gas cylinders for tcpO ₂ /tcpCO ₂ calibration. (meets international shipping requirements)	6 cylinders	15205A, 15206A, 15209A, 15209AH, 15210A, 15210AH, 15210B, M1018A, M1918A, M1918B, M3015A
15210-64020	Gas Cal 2 gas cylinders for tcpO ₂ /tcpCO ₂ calibration. (meets international shipping requirements)	6 cylinders	15205A, 15206A, 15209A, 15209AH, 15210A, 15210B, 15210AH, M1018A, Microstream CO ₂ Ext. M3015A
M2205A	Calibration tubes used to connect 15210B calibration unit to M1918A or M1918B transducer during calibration.	5 tubes	15210B, M1018A, M1918A M1918B

ECG Cables and Leads

3-Lead IntelliVue ECG Patient Cable System

- Replaces CMS cable system as well as individual lead sets
- Lightweight and comfortable for tiny patients

Product #	Description	Packaging	Use With:
M1669A CE	3-Lead ECG Trunk, AAMI/IEC (12 pin connector) Length: 9' (2.7m)	1 cable	Multi-Measurement Modules: M3000A, M3001A, M3002A, ECG Modules: M1001B, M1002B, MP5 Patient Monitor: M8105A, MP2 Patient Monitor: M8102A, SureSigns Monitors VM4: 863063, VM6: 863064, 863065, VM8: 863066, 863068
M1622A	3-Lead Miniclip lead set, AAMI Length: 1.5 feet, grabbers, non-shielded (purple)	1 lead set	Use with M1669A Trunk Cable
M1624A	3-Lead Miniclip lead set, AAMI Length: 2.3 feet, grabbers, non-shielded (purple)	1 lead set	Use with M1669A Trunk Cable
M1626A CE	3-Lead Miniclip lead set, IEC Length: 0.7m, grabbers, non-shielded (purple)	1 lead set	Use with M1669A Trunk Cable
989803170171 CE	CBL OR 3-Lead ECG Trunk Cable, AAMI/IEC, Cable length = 2.7m (9')	1 cable	Compatible Bedside Monitors M1001B, M1002B, M3000A, M3001A, M3002A, M8105A (MP5), M8102A (MP2), 863063 (VM4), 863064 (VM6), 863065 (VM6), 863066 (VM8), 863068 (VM6), 863077 (VM3), Compatible Defibrillators M3535A (MRx monitor/defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL)



M1669A



M1622A



M1624A



M1626A

5-Lead IntelliVue ECG Patient Cable System

- Replaces CMS cable system as well as individual lead sets
- Lightweight and comfortable for tiny patients

Product #	Description	Packaging	Use With:
M1668A CE	5-Lead ECG Trunk AAMI/IEC (12 pin connector) Length: 9 feet (2.7m)	1 cable	ECG Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5 Patient Monitor: M8105A, MP2 Patient Monitor: M8102A, SureSigns Monitors: VM6: 863064, 863065, VM8: 863066, 863068, 860338, 860335
M1647A	5-Lead Miniclip lead set, AAMI Length: 2.3 feet/4.2 feet, grabbers, non-shielded (purple)	1 lead set	Use with M1668A Trunk Cable
M1648A CE	5-Lead Miniclip lead set, IEC Length: 0.7m/1.3m, grabbers, non-shielded (purple)	1 lead set	Use with M1668A Trunk Cable
989803170181 CE	CBL OR 5-Lead ECG Trunk Cable, AAMI/IEC, Cable length = 2.7m (9ft)	1 cable	Compatible Bedside Monitors M1001B, M1002B, M3000A, M3001A, M3002A, M8105A (MP5), M8102A (MP2), 863064 (VM6), 863065 (VM6), 863066 (VM8), 863068 (VM6), Compatible Defibrillators M3534A (MRx monitor/defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL)



M1668A



M1647A



M1648A

3-Lead CMS ECG Patient Cable System - AAMI

Product #	Description	Packaging	Use With:
M1500A	3-Lead ECG Trunk, AAMI (12 pin connector) Length: 9'	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, C3 Monitors: 862474, 862478, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077
M1540C	3-Lead ECG Trunk, AAMI (12 pin connector) Length: 3'	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A
M1608A	3-Lead Miniclip lead set, AAMI Length: 1.5', non-shielded (purple)	1 lead set (3 leads)	Use with M1500A or M1540C Trunk Cable
M1609A	3-Lead Miniclip lead set, AAMI Length: 2.3', non-shielded (purple)	1 lead set (3 leads)	Use with M1500A or M1540C Trunk Cable

3-Lead CMS ECG Patient Cable System - IEC

Product #	Description	Packaging	Use With:
M1510A CE	3-Lead ECG Trunk, IEC (12 pin connector) Length: 2.7m	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, C3 Monitors: 862474, 862478, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077
M1619A CE	3-Lead Miniclip lead set, IEC Length: 0.7m, non-shielded (purple)	1 lead set (3 leads)	Use with M1510A Trunk Cable

5-Lead CMS ECG Patient Cable System - AAMI

Product #	Description	Packaging	Use With:
M1520A	5-Lead ECG Trunk Cable, AAMI (12 pin connector) Length: 9'	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, C3 Monitors: 862474, 862478, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077
M1560C	5-Lead ECG Trunk Cable, AAMI (12 pin connector) Length: 2.9'	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A
M1629A	5-Lead Miniclip lead set, AAMI Length: 2.3'/4.2', grabbers, non-shielded (purple)	1 lead set (5 leads)	Use with M1520A or M1560C Trunk Cable

5-Lead CMS ECG Patient Cable System - IEC

Product #	Description	Packaging	Use With:
M1530A CE	5-Lead ECG Trunk Cable, IEC (12 pin connector) Length: 2.7m	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, C3 Monitors: 862474, 862478, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077
M1570C CE	5-Lead ECG Trunk Cable, IEC (12 pin connector) Length: 0.9m	1 cable	ECG Modules: M1001A, M1001B, M1002A, M1002B, Legacy HP Monitors: 78352A, 78352C, 78354C, 78834C, IntelliVue Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A
M1639A CE	5-Lead Miniclip lead set, IEC Length: 0.7m/1.3m, grabbers, non-shielded (purple)	1 lead set (5 leads)	Use with M1530A or M1570C Trunk Cable

Radiolucent ECG Lead Sets – AAMI/IEC Color Coding

Product #	Description	Packaging	Use With:
M1649A CE	Radiotranslucent unshielded 5-lead set with safety connectors. Grabber termination, 5 individual wires, IEC/AAMI color coding Length: 36" (914mm)	5 individual leads	MMM X2: M3002A; MP2: M8102A; MP5: M8105A, 860335, 860338



M1649A

ECG Electrodes


Disposable Lead Wire Electrodes: AAMI Color Coding

Metallic - AAMI


Product #	Description	Packaging	Use With:
13952A	Preattached metallic (18") leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 1-1/4" (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors
13953A CE	Preattached metallic (18") leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 7/8" (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors

Radiolucent - AAMI

Product #	Description	Packaging	Use With:
13952B CE	Preattached radiolucent (18") leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 1-1/4" (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors
13953B	Preattached radiolucent (18") leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 7/8" (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors



13952A



13952B

Disposable Lead Wire Electrodes: IEC Color Coding

Metallic - IEC

Product #	Description	Packaging	Use With:
13952D CE	Preattached metallic (46cm) leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 33mm, (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors
13953D CE	Preattached metallic (46cm) leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 22mm, (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors

Radiolucent - IEC

Product #	Description	Packaging	Use With:
13952E CE	Preattached radiolucent (46cm) leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 33mm, (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors
13953E CE	Preattached radiolucent (46cm) leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 22mm, (socket size: 1.5mm), shelf life: 6 months	3 electrodes/pack 300/case	Pediatric/Neonatal Monitors



13953D

Disposable Snap Connector Electrodes

Product #	Description	Packaging	Use With:
13950B CE	Disposable cloth electrode with KCl foam gel retaining ring, snap. (Individually wrapped) Diameter: 1.25" (32mm), shelf life: 6 months	30 electrodes/box 10 boxes/case	Pediatric/Neonatal Monitors
13951C CE	Disposable, repositionable solid gel electrode, snap, shelf life: 6 months. (3 electrodes mounted on plastic card) Dimensions: 7/8" x 1-1/4" (22mm x 33mm)	30 electrodes/pack 300/case	Pediatric/Neonatal Monitors
13955C CE	Disposable, repositionable solid gel electrode, snap, shelf life: 6 months. (3 electrodes mounted on plastic card) Dimensions: 7/8" x 7/8" (22mm x 22mm)	30 electrodes/pack 300/case	Pediatric/Neonatal Monitors



13950B



13951C



13955C

Invasive Blood Pressure

TRANSPAC® 4 Neonatal Disposable Pressure Transducer Monitoring Kits

Product #	Description	Packaging	Use With:
989803179841 CE 0123	Single Neonatal DPT TP4 24" (61cm), Two 3-Way Lipid Resistant Polysulfone Stopcocks, 30ml/hr Flush Device	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A
989803179851 CE 0123	Single Neonatal DPT TP4 12" (31cm), Three 3-Way Lipid Resistant Polysulfone Stopcocks, 30ml/hr Flush Device	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A
989803179881 CE 0123	Single Neonatal DPT TP4 18" (46cm), Three 3-way Lipid Resistant Polysulfone Stopcocks, 30ml/hr Flush Device, One in-line Luer Sampling Port and 10ml Reservoir	20 kits	M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A



989803179841



989803179881

TRANSPAC® 4 Reusable Transducer Cables

Product #	Description	Packaging	Use With:
989803177921 CE0123	Transpac 4 Cable for Philips Monitors 10' (3m)	1 cable	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179941 CE0123	Transpac 4 Cable for Philips Monitors 15' (4.5m)	1 cable	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A
989803179951 CE0123	Transpac 4 Cable for Philips Monitors 15' (4.5m), bulk package	20 cables	TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A

TRANSPAC® 4 Disposable Transducer Kit Accessories

Product #	Description	Packaging	Use With:
989803177941 CE 0123	Transducer Holder for IV Pole, Reusable	5 holders	TRANSPAC 4 Disposable Transducer Kits
989803179911 CE 0123	Single Transducer Mount, Reusable	30 mounts	TRANSPAC 4 Disposable Transducer Kits
989803180861 CE 0123	Transducer Simulator, Reusable	1 simulator	TRANSPAC 4 Disposable Transducer Kits and Cables with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A



989803177941



989803179911

Non-Invasive Blood Pressure

Reusable NIBP Cuffs

Reusable NIBP Cuffs

Product #	Description	Packaging	Use With:
M1571A CE	Multi-Patient Comfort Cuff, 1 Hose, Infant Bladder width: 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10cm to 15cm	1 cuff	NIBP Modules: M1008B, Multi-Measurement Server: M3000A, M3001A, M3002A; MP5: M8105A, MP2: M8102A, C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 866077, M8105AS, 860335, 860338, 863275, 863276, 863278
M4552B CE	Easy Care Cuff, 1 Hose, Infant Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm Cuff length = 8.7" (22cm)	1 cuff	NIBP Modules: M1008B, Multi-Measurement Server: M3000A, M3001A, M3002A; MP5: M8105A, MP2: M8102A, C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 866077, M8105AS, 863275, 863276, 863278
M4552B5 CE	Easy Care Cuff Kit, 1 Hose, Infant Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm Cuff length = 8.7" (22cm)	5 cuffs	NIBP Modules: M1008B, Multi-Measurement Server: M3000A, M3001A, M3002A; MP5: M8105A, MP2: M8102A, C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 866077, M8105AS



M1571A



M4552B

Disposable NIBP Cuffs

Soft Single Patient Neonatal/Infant Disposable Cuffs

Product #	Description	Packaging	Use With:
M1866S CE	Neonatal Soft Single-Patient Cuff Size #1 Disposable, single hose	40 cuffs per box	NIBP Modules, M1008B, Multi-Measurement Modules, M3000A, M3001A, M3002A, MP5, M8105A, MP2, M8102A, SureSigns Monitors, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
M1868S CE	Neonatal Soft Single-Patient Cuff Size #2 Disposable, single hose	40 cuffs per box	NIBP Modules, M1008B, Multi-Measurement Modules, M3000A, M3001A, M3002A, MP5, M8105A, MP2, M8102A, SureSigns Monitors, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
M1870S CE	Neonatal Soft Single-Patient Cuff Size #3 Disposable, single hose	40 cuffs per box	NIBP Modules, M1008B, Multi-Measurement Modules, M3000A, M3001A, M3002A, MP5, M8105A, MP2, M8102A, SureSigns Monitors, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
M1872S CE	Neonatal Soft Single-Patient Cuff Size #4 Disposable, single hose	40 cuffs per box	NIBP Modules, M1008B, Multi-Measurement Modules, M3000A, M3001A, M3002A, MP5, M8105A, MP2, M8102A, SureSigns Monitors, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278

Soft Single Patient Neonatal/Infant Disposable Cuffs (continued)

Product #	Description	Packaging	Use With:
M1873S CE	Infant Soft Single-Patient Cuff Size #5 Disposable, single hose	40 cuffs per box	NIBP Modules, M1008B, Multi-Measurement Modules, M3000A, M3001A, M3002A, MP5, M8105A, MP2, M8102A, SureSigns Monitors, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
989803167551 CE	Neonatal Soft Single-Patient NIBP Cuffs Disposable, single hose (consists of 5 x M1866S, 10 x M1868S, 20 x M1870S, 10 x M1872S & 5 x M1873S)	50 cuffs per box	NIBP Modules, M1008B, Multi-Measurement Modules, M3000A, M3001A, M3002A, MP5, M8105A, MP2, M8102A, SureSigns Monitors, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS



M1866S, M1868S, M1870S, M1872S, M1873S

Single Patient Neonatal/Infant Disposable Cuffs

Product #	Description	Packaging	Use With:
M1866B CE	Neonatal Single Patient Cuff, Size 1 Circumference: 3.1 – 5.7cm, 1 Hose. New connector configuration, replaces the M1866A.	40 cuffs	NIBP Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
M1868B CE	Neonatal Single Patient Cuff, Size 2 Circumference: 4.3 – 8.0cm, 1 Hose. New connector configuration, replaces the M1868A.	40 cuffs	NIBP Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
M1870B CE	Neonatal Single Patient Cuff, Size 3 Circumference: 5.8 – 10.9cm, 1 Hose. New connector configuration, replaces the M1870A.	40 cuffs	NIBP Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
M1872B CE	Neonatal Single Patient Cuff, Size 4 Circumference: 7.1 – 13.1cm, 1 Hose. New connector configuration, replaces the M1872A.	40 cuffs	NIBP Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278

<p>M1873B CE</p>	<p>Infant Single Patient Cuff, Size 5 Circumference: 10.0 – 15.0cm, 1 Hose New connector configuration.</p>	<p>40 cuffs</p>	<p>NIBP Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278</p>
<p>989803167541 CE</p>	<p>Neonatal Single Patient NIBP Kit contains 50 Cuffs. Incorporates new connector configuration and replaces M1820-60020.</p>	<p>5 M1866B Size 1, 10 M1868B Size 2, 20 M1870B Size 3, 10 M1872B Size 4, 5 M1873B Size 5</p>	<p>NIBP Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS</p>



M1866B, M1868B, M1870B, M1872B, M1873B

Single Patient Disposable NIBP Cuffs

Product #	Description	Packaging	Use With:
REVISED M1820-60020	Neonatal Single Patient NIBP Cuff Kit, Sizes 1 – 4	10 M1866A Size 1, 10 M1868A Size 2, 20 M1870A Size 3, 10 M1872A Size 4.	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5 M8105A, MP2 M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors: 78352A, 78352C, 78834C
REVISED M1866A	Neonatal Single Patient NIBP Cuff, Size 1	20 cuffs	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
REVISED M1868A	Neonatal Single Patient NIBP Cuff, Size 2	20 cuffs	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
REVISED M1870A	Neonatal Single Patient NIBP Cuff, Size 3	20 cuffs	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278



M1866A



M1868A



M1870A

Single Patient Disposable NIBP Cuffs *(continued)*

Product #	Description	Packaging	Use With:
M1872A	Neonatal Single Patient NIBP Cuff, Size 4	20 cuffs	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP5: M8105A, MP2: M8102A, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
M1874A CE	Single Patient NIBP Cuff, Infant	10 cuffs	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5: M8105A, MP2: M8102A, C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082, M8105AS, 860335, 860338, 863275, 863276, 863278
M4572B CE	Soft Single Patient Disposable NIBP Cuff, Infant	10 cuffs	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5: M8105A, MP2: M8102A, C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M3535A, M3536A, M8102A, M8105A, M8105AS, M8105AT, 860335, 860338, 863275, 863276, 863278



M1872A



M1874A



M4572B

Blood Pressure Air Hoses

Product #	Description	Packaging	Use With:
M1596C CE	Neonatal Blood Pressure Air Hose. Length = 4.92' (1.5m). Connects to Philips Neonatal Single Patient Cuffs only and incorporates new connector configuration. This cable is used for close patient proximity and replaces the M1596B.	1 air hose	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5 M8105A, MP2: M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
M1597C CE	Neonatal Blood Pressure Air Hose. Length = 9.84' (3.0m). Connects to Philips Neonatal Single Patient Cuffs only and incorporates new connector configuration. This cable is used for classic bedside monitoring and replaces the M1597B.	1 air hose	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5 M8105A, MP2: M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
989803167561 CE	Neonatal NIBP Air Hose for use with the Philips Vital Signs Monitor. Incorporates new connector configuration and replaces the 989803136941.	1 air hose	VS1 Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062
989803167511 CE	Connector upgrade kit for M1597B. Upgrade existing NIBP airhoses with an inner diameter of 2mm. Contains 20 female connectors with 2mm barb and instructions.	20 connectors	Upgrades M1597B to M1597C
989803167521 CE	Connector upgrade kit for M1596B. Upgrade existing NIBP airhoses with an inner diameter of 3mm. Contains 20 female connectors with 3mm barb and instructions.	20 connectors	Upgrades M1596B to M1596C



M1596C



M1597C



989803167561



989803167511



989803167521

Blood Pressure Air Hoses (continued)

Product #	Description	Packaging	Use With:
989803167531 CE	Connector upgrade kit for 989803136941. Upgrade existing NIBP airhoses with an inner diameter of 4mm. Contains 20 female connectors with 4mm barb and instructions.	20 connectors	Upgrades 989803136941 Neonatal Hose to 989803167561
989803168821 CE	Connector upgrade kit for Philips Neonatal Cuffs. Contains 20 male connectors and instructions.	20 connectors	Upgrades M1866A, M1868A, M1870A, and M1872A to work with M1596C and M1597C hoses
989803136941	Neonatal Blood Pressure Air Hose (not available in Europe) Length: 11.71' (3.57m)	1 hose	VS1 Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062
M1596B	Neonatal Blood Pressure Air Hose (not available in Europe) length: 4.92' (1.5m)	1 hose	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5 M8105A, MP2 M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
M1597B	Neonatal Blood Pressure Air Hose (not available in Europe) length: 9.84' (3.0m)	1 hose	NIBP Modules: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP5 M8105A, MP2 M8102A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278



M1597B

Pulse Oximetry

Single Patient Disposable SpO₂ Sensors

Philips Disposable SpO₂ Sensors

Product #	Description	Packaging	Use With:
M1132A CE	Infant Disposable Wrap Sensor Infant finger. Patient size: >6.6 – 22lbs. (>3 – 10kg)	20 sensors	Use Adapter Cable M1900B for SpO ₂ Module: M1020A. Use Adapter Cable M1943A or M1943AL for SpO ₂ Module: M1020B Opt A01;
M1133A CE	Neonatal/Infant/Adult Disposable Wrap Sensor Neonatal foot/hand. Patient size: <6.6lbs. (<3kg) Infant big toe/thumb. Patient size: 22 – 44lbs. (10 – 20kg)	20 sensors	Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01; MP5 Monitor: M8105A, MP2: M8102A; SureSigns Monitors: VM1: 863264, 863265, 863266, 863278
M1134A CE	Adhesive-Free Neonatal/Infant/Adult Disposable Wrap Sensor. Neonatal foot/hand. Patient size: <6.6lbs. (<3kg) Infant big toe/thumb. Patient size: 22 – 44lbs. (10 – 20kg)	20 sensors	VM3/VM4/VM6/VM8: 863077, 863063, 863064, 863065, 863066, 863068, VS2: 863080, 863081, 863082; VS3: 863069, 863071, 863072, 863073, 863074, M8105AS, 860335, 860338



M1132A



M1133A



M1134A

Adhesive SpO₂ Sensors from Philips

Product #	Description	Packaging	Use With:
M1901B CE	Neonatal Disposable SpO ₂ Sensor Neonatal foot. Patient size: <7 lbs. (<3kg) (not available in N. America or Japan)	24 sensors	Use Adapter Cable M1943A or M1943AL for SpO ₂ Modules: M1020B Opt. A01; Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01; M1020B-A04 MP5 Monitor: M8105A, MP2: M8102A; M8105AS Suresigns Monitors: 863264, 863265, 863266, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863276, 863278
M1902B CE	Infant Disposable SpO ₂ Sensor Infant digit. Patient size: 7lbs. to 44lbs. (3 – 20kg) (not available in N. America or Japan)	24 sensors	Use Adapter Cable M1943NL for SpO ₂ Modules: M1020B Opt. A02; Multi-Measurement Modules: M3001A Opt. A02 or Opt. A04; M3002A Opt. A02 or A04; MP5 Monitor: M8105A Opt. SP4; MP2 Monitor: M8102A Opt. SP4, M8105AS, 863276, 863278



M1901B



M1902B

Reusable SpO₂ Sensors

Philips Reusable SpO₂ Sensors

Product #	Description	Packaging	Use With:
REVISED M1193A CE	Reusable Neonatal Wrap Sensor Length = 1.5m Neonatal foot/hand. Patient size: 2.2lbs. to 8.8lbs. (1-4kg) 8-pin connector	1 sensor	Use Adapter Cable M1940A for SpO ₂ Modules: M1020A Use Extension Cable M1941A for SpO ₂ Modules: M1020B Opt A01; Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01; MP5 Monitor: M8105A, MP2: M8102A; Suresigns Monitors: VM1: 863264, 863265, 863266, , 863276, 863278, VM3/VM4/VM6/VM8: 863077, 863063, 863064, 863065, 863066, 863068, VS2: 863080, 863081, 863082; VS3: 863071, 863072, 863073, 863074, M8105AS
REVISED M1193T CE	Reusable Neonatal Wrap Sensor Length = 1.5m Neonatal foot/hand. Patient size: 2.2lbs. to 8.8lbs. (1 – 4kg) 9-pin (D-sub) connector.	1 sensor	Use Adapter Cable M1900B for SpO ₂ Modules: M1020A Use Adapter Cable M1943A or M1943AL for SpO ₂ Modules: M1020B Opt A01; Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01; MP5 Monitor: M8105A, MP2: M8102A; Suresigns Monitors: VM1: 863264, 863265, 863266, 863276, 863278, VM3/VM4/VM6/VM8: 863077, 863063, 863064, 863065, 863066, 863068, VS2: 863080, 863081, 863082; VS3: 863071, 863072, 863073, 863074, M8105AS



M1193A



M1193T

Philips Reusable SpO₂ Sensors *(continued)*

Product #	Description	Packaging	Use With:
<p>REVISED</p> <p>M1195A</p> <p>CE</p>	<p>Reusable Infant Finger Glove Sensor</p> <p>Length = 1.5m</p> <p>Infant finger. Patient size: 8.8lbs. to 33lbs. (4 – 15kg)</p> <p>8-pin connector</p>	<p>1 sensor</p>	<p>Use Adapter Cable M1940A for SpO₂ Modules: M1020A;</p> <p>Use Extension Cable M1941A for SpO₂ Modules: M1020B Opt A01;</p> <p>Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01;</p> <p>MP5 Monitor: M8105A, MP2: M8102A;</p> <p>SureSigns Monitors: VM1: 863264, 863265, 863266, , 863276, 863278, VM3/VM4/VM6/ VM8: 863077, 863063, 863064, 863065, 863066, 863068, VS2: 863080, 863081, 863082; VS3: 863071, 863072, 863073, 863074, M8105AS</p>



M1195A

SpO₂ Adapter Cables

Product #	Description	Packaging	Use With:
<p>REVISED</p> <p>989803148221</p> <p>CE</p>	<p>Masimo® MP10 Adapter Cable (LNCS sensors)</p> <p>Length: 10' (3m)</p>	1 cable	<p>SpO₂ Modules: M1020B Opt A01, M1020B Opt A03;</p> <p>Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01, M3001A Opt. A03, M3002A Opt. A03;</p> <p>MP5 Monitor: M8105A, MP2: M8102A</p> <p>SureSigns Monitors: 863264, 863265, 863266, 863077, 863063, 863064, 863065, 863066, 863068, 863080, 863081, 863082; 863071, 863072, 863073, 863074, M8105AS</p>
<p>M1020-61100</p> <p>CE</p>	<p>Masimo® MP12 Adapter Cable (LNOP sensors)</p> <p>Length: 11.8' (3.6m)</p>	1 cable	<p>SpO₂ Modules: M1020B Opt A01, M1020B Opt A03;</p> <p>Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01, M3001A Opt. A03, M3002A Opt. A03;</p> <p>MP5 Monitor: M8105A, MP2: M8102A</p> <p>SureSigns Monitors: 863264, 863265, 863266, 863077, 863063, 863064, 863065, 863066, 863068, 863080, 863081, 863082; 863071, 863072, 863073, 863074, 863276, 863278</p>
<p>REVISED</p> <p>M1020-61102</p> <p>CE</p>	<p>Masimo LNOP Adapter Cable (12-pin)</p> <p>Adapts LNOP sensors to Viridia, CMS monitors</p> <p>Length: 11.8' (3.6m)</p>	1 cable	SpO ₂ Module: M1020A
<p>M1900B</p> <p>CE</p>	<p>SpO₂ Adapter Cable (12-pin)</p> <p>Adapts 9-pin sensors to 12-pin sockets</p> <p>Length: 9.8' (3m)</p>	1 cable	SpO ₂ Module: M1020A



989803148221



M1020-61100



M1900B

SpO₂ Adapter Cables *(continued)*

Product #	Description	Packaging	Use With:
M1940A CE	SpO ₂ Adapter Cable (12-pin) Adapts 8-pin sensors to 12-pin sockets Length: 6.6' (2m)	1 cable	SpO ₂ Modules: M1020A
M1941A CE	SpO ₂ Extension Cable (8-pin) Adds additional length to 8-pin sensors to 8-pin sockets Length: 6.6' (2m)	1 cable	SpO ₂ Modules: M1020B Opt A01; Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01; MP5 Monitor: M8105A, MP2: M8102A; M8105AS SureSigns Monitors: VM1: 863264, 863265, 863266, VM3/VM4/VM6/VM8: 863077, 863063, 863064, 863065, 863066, 863068, VS2: 863080, 863081, 863082; VS3: 863071, 863072, 863073, 863074
M1943A CE	SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 3.3' (1.1m)	1 cable	SpO ₂ Modules: M1020B Opt A01; Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, M3002A Opt. A01; MP5 Monitor: M8105A, MP2: M8102A; M8105AS, 860335, 860338
M1943AL CE	SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 9.8' (3m)	1 cable	SureSigns Monitors: VM1: 863264, 863265, 863266, VM3/VM4/VM6/VM8: 863077, 863063, 863064, 863065, 863066, 863068, VS2: 863080, 863081, 863082; VS3: 863071, 863072, 863073, 863074, M8105AS
REVISED M1943NL CE	OxiMax® SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to IntelliVue A02/A04 sockets Length: 9.8' (3m)	1 cable	M1020B-A02, M3001A-A02, M3001A-A04, M3002A-A04, M8102A-SP4, M8105A-SP4, M8105AT-SP4, M8105AS, M1020B-A04



M1940A



M1941A



M1943A/M1943AL



M1943NL

Temperature

Disposable Skin Surface Temperature Probe

Product #	Description	Packaging	Use With:
21091A CE	Disposable sterile skin surface temperature probe Length: 30" (76cm) Shelf life: at least 12 months	1 box = 20 probes	Temperature Module: M1029A Multi-Measurement Modules: M3000A, M3001A, M3002A SureSigns: 863064, 863065, 863066, 863068, M3015A
21082A CE	Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 10' (3m)	1 cable	MP2 Monitor: M8102A, MP5 Monitor: M8105A, M3015A
21082B CE	Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 5' (1.5m)	1 cable	M3015A
989803162591 CE	Disposable temperature probe adapter cable for Series 400 monitors Length: 10' (3m)	1 cable	Series 400 Temperature Monitors with 1/4-inch phone plug
989803162601 CE	Disposable temperature probe adapter cable for Series 400 monitors Length: 5' (1.5m)	1 cable	Series 400 Temperature Monitors with 1/4-inch phone plug



21091A



21082A/B



989803162591



989803162601

Reusable Skin Surface Temperature Probe

Product #	Description	Packaging	Use With:
21078A CE	Reusable skin surface temperature probe, autoclaveable (2-pin connector) Diameter: 0.37" (9.5mm) Length: 10' (305cm)	1 probe	No cable require for: Temperature Module: M1029A Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A SureSigns: 863064, 863065, 863066, 863068, 860335, 860338 MP2 Monitor: M8102A, MP5 Monitor: M8105A
989803162641 CE	Reusable skin surface temperature probe, autoclaveable (1/4-inch phone plug connector) Diameter: 0.37" (9.5mm) Length: 10' (305cm)	1 probe	No cable required for: Series 400 Temperature Monitors
989803163271 CE	Temperature probe adapter cable. Required to connect 21078A (2-pin) temperature probe to Series 400 temperature monitors with 1/4" phone plug	1 probe	Series 400 Temperature Monitors with 1/4 inch phone plug
989803163281 CE	Temperature probe adapter cable. Required to connect 989803162641 (1/4" phone plug) temperature probe to Philips and legacy HP/Agilent temperature monitors with 2-pin input connector.	1 probe	Philips and legacy HP/Agilent monitors with 2-pin input connector



21078A



989803162641



989803163271




989803163281

Index

13950B	.46	989803137641	.13	989803178021	.32	M1194A	.27
13951C	.46	989803137651	.12	989803178031	.32	M1195A	.62
13952A	.44	989803137791	.12	989803178041	.32	M1196A	.27
13952B	.44	989803137801	.12	989803178051	.32	M1196S	.27
13952D	.45	989803137811	.12	989803178071	.32	M1196T	.27
13952E	.45	989803139771	.12	989803179121	.32	M1333A	.15
13953A	.44	989803143401	.9	989803179841	.47	M1334A	.15
13953B	.44	989803143411	.16	989803179851	.47	M1363A	.18
13953D	.45	989803143931	.15	989803179881	.47	M1500A	.41
13953E	.45	989803144471	.33	989803179911	.49	M1510A	.41
13955C	.46	989803144531	.34	989803179941	.48	M1520A	.42
14328A	.9	989803148221	.63	989803179951	.48	M1530A	.42
1500-0628	.10	989803151601	.25	989803180861	.49	M1540C	.41
1500-0643	.10	989803151951	.37	989803182291	.24	M1560C	.42
15155-43101	.9	989803151961	.37	989803182301	.24	M1562A	.10
15209-60010	.37	989803159571	.31	989803182311	.24	M1562B	.11
15209-60020	.37	989803159581	.31	989803182321	.24	M1569A	.9
15210-60010	.38	989803159581	.31	989803183331	.22	M1570C	.42
15210-60020	.38	989803160241	.31	989803183341	.22	M1571A	.50
15210-64010	.38	989803160251	.31	989803183351	.22	M1573A	.19
15210-64020	.38	989803160261	.31	989803183361	.22	M1574A	.19
		989803160261	.31	989803183371	.22	M1575A	.19
21078A	.66	989803160271	.31			M1576A	.19
21082A	.65	989803160371	.37	CPJ840J5	.15	M1578A	.21
21082B	.65	989803162591	.65	CPJ84022	.15	M1596B	.58
21091A	.65	989803162601	.65	CPJ84026	.15	M1596C	.57
		989803162641	.66	CPJ84046	.15	M1597B	.58
40483A	.18	989803163271	.66	CPJ84447	.15	M1597C	.57
40483B	.18	989803163281	.66			M1598B	.25
40493D	.18	989803167511	.57	M1020-61100	.63	M1599B	.25
40493E	.18	989803167521	.57	M1020-61102	.63	M1608A	.41
		989803167531	.58	M1131A	.28	M1609A	.41
862214	.17	989803167541	.54	M1132A	.59	M1619A	.41
862215	.17	989803167551	.52	M1133A	.28	M1622A	.39
862217	.17	989803167561	.57	M1133A	.59	M1624A	.39
862218	.17	989803168821	.58	M1134A	.28	M1626A	.39
864288	.21	989803170171	.39	M1134A	.59	M1627A	.27
864290	.21	989803170181	.40	M1191B	.26	M1629A	.42
864291	.21	989803177921	.48	M1191BL	.26	M1639A	.42
		989803177941	.49	M1191T	.26	M1647A	.40
9270-0485	.14	989803177971	.31	M1192A	.26	M1648A	.40
9270-0630-050	.14	989803177981	.31	M1192T	.26	M1649A	.43
989803136941	.58	989803178001	.31	M1193A	.61	M1668A	.40
989803137631	.13	989803178011	.31	M1193T	.61	M1669A	.39

M1820-60020	.55	M2524A	.32
M1866A	.55	M2533A	.35
M1866B	.53	M2536A	.35
M1866S	.51	M2738A	.16
M1868A	.55	M2741A	.33
M1868B	.53	M2745A	.33
M1868S	.51	M2746A	.33
M1870A	.55	M2751A	.33
M1870B	.53	M2757A	.33
M1870S	.51	M2761A	.34
M1872A	.56	M2773A	.34
M1872B	.53	M2776A	.34
M1872S	.51	M2777A	.34
M1873B	.54	M2781A	.36
M1873S	.52	M2782A	.36
M1874A	.56	M2783A	.36
M1900B	.29	M2785A	.36
M1900B	.63	M2786A	.36
M1901B	.60	M4503A	.16
M1902B	.60	M4552B	.50
M1904B	.28	M4552B5	.50
M1910A	.14	M4554B/ M4554B5	.20
M1911A	.14	M4555B/ M4555B5	.20
M1913A	.14	M4556B/ M4556B5	.20
M1918A	.37	M4557B/ M4557B5	.20
M1918B	.37	M4558B/ M4558B5	.20
M1920A	.32	M4559B/ M4559B5	.20
M1921A	.32	M4572B	.56
M1923A	.32	M4574B	.23
M1940A	.29	M4575B	.23
M1940A	.64	M4576B	.23
M1941A	.29	M4577B	.23
M1941A	.64	M4578B	.23
M1943A	.29	M4579B	.23
M1943A	.64	M4601A	.11
M1943AL	.29	M4602A	.11
M1943AL	.64	M4603A	.11
M1943NL	.64	M4681A	.32
M2205A	.38	M4687A	.32
M2208A	.9	M4691A	.32
M2501A	.35		
M2513A	.35		
M2516A	.35		
M2520A	.32		



© 2013 Koninklijke Philips N.V. All rights are reserved.
Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands
4522 962 95511 * JUN 2013

How to reach us
www.philips.com/healthcare
healthcare@philips.com

Product information
www.philips.com/healthcarestore