



A New Design
Reducing the Risk of AV Block

1 Reducing the Risk of AV Block

Elongatable waist designed for a variety of VSD

Soft woven mesh provides high conformability to septal defects

2 Low Profile



Partnered with 4-7 Fr. delivery sheath


Slim cable to minimize unwanted damage

3 Procedural Flexibility

Double-sided screw for retrograde or antegrade approach

Ordering Information

Catalogue Number	D Disc Diameter (mm)	D2 Waist Diameter LV Side (mm)	D1 Waist Diameter RV Side (mm)	L Waist Length (mm)	Recommended Delivery Sheath
	LT-MFO-5-3	10	5	3	4F-5F
	LT-MFO-6-4	10	6	4	
	LT-MFO-7-5	12	7	5	
	LT-MFO-8-6	12	8	6	
 with membrane	LT-MFO-9-7	14	9	7	6F
	LT-MFO-10-8	14	10	8	
	LT-MFO-12-10	16	12	10	7F
	LT-MFO-14-12	18	14	12	



Lifetech Scientific (Shenzhen) Co., Ltd.

Add: 8F, LifeTech Scientific Building, No.22, Keji 12th Road South, High-tech Industrial Park, Yuehai Subdistrict, Nanshan District, Shenzhen 518063, P. R. China.

Tel: +86-755-86026250 Fax: 86-755-86026251

Email: marketing@lifetechmed.com

Web: www.lifetechmed.com

EU authorized representative:
LifeTech Scientific (Europe) Coöperatief U.A
Oliemolenstraat 60, 6416CB Heerlen, the Netherlands

2021-EN-2.0