



www.vathin.com
marketing@vathin.com

1st floor, Building 12 , No.9, West Chuanqi Road, Jiu Hua
Economic and Technological Development Zone, Xiangtan
411100, China



Single-use Flexible Video Bronchoscope H-SteriScope™ Series

Diverse Options, Unified Excellence


Key Benefits

H-SteriScope Series offers seven different models, accommodating to various diagnostic and therapeutic needs. Zero and Ultra Slim models are available for pediatric procedures and Extra Large is designed with instrument-intensive applications.

This innovative technology enhances the scope's stability, durability and flexibility, ensuring precise maneuvering and safe operation.

Designed for optimal portability and ease of use, the device supports diverse clinical applications, enabling procedures in environments such as emergency rooms, intensive care units, and patient bedsides.

By eliminating reprocessing, single-use bronchoscopes reduce cross-infection risks and avoid wear and tear, enhancing patient safety and cost efficiency.



Digital Video Monitor	
Model	DVM-B1
Video Record	Duration more than 6 hours
Playback	Picture & Video
Electrical Power Requirement	100V-240V AC; 47-63HZ; 1.62-0.72A
Dimension WxHxT	310mmx231mmx39mm
Weight	2.14kg
Storage Capacity	64G



7 Models, 1 Solution



Model	ID	OD
Zero	-	2.2mm
Ultra Slim	1.2mm	2.6mm
Slim	1.2mm	3.2mm
IP Scope	2.0mm	4.0mm
Normal	2.2mm	4.9mm
Large	2.8mm	5.8mm
Extra	3.2mm	6.2mm

Specifications

Name	Zero	Ultra Slim	Slim	IP Scope	Normal	Large	Extra
Model	BCV1-02	Ultra Slim	BCV1-C2	BCV1-KS3-L	BCV1-M2	BCV1-S2	BCV1-W2
Bending Section	U210°/D210°	U210°/D210°	U210°/D210°	U210°/D210°	U210°/D210°	U210°/D210°	U210°/D210°
Insertion Tube OD	2.2mm	2.6mm	3.2mm	4.0mm	4.9mm	5.8mm	6.2mm
Working channel	-	1.2mm	1.2mm	2.0mm	2.2m	2.8mm	3.2mm
Working Length	600mm			620mm	600mm		
Field of view	110°±5%						
Depth of view	3~100mm						
Illumination method	LED						
Suction Connector	Ø 7± 1						
Method of Sterilization	ETO						